

Test Definition: BCRFX

BCR/ABL1 Qualitative Diagnostic Assay with
Reflex to BCR/ABL1 p190 Quantitative Assay or
BCR/ABL1 p210 Quantitative Assay, Varies

Overview

Useful For

Diagnostic workup of patients with high probability of *BCR::ABL1*-positive hematopoietic neoplasms, predominantly chronic myeloid/myelogenous leukemia and acute lymphoblastic leukemia

Reflex Tests

| Test Id | Reporting Name | Available Separately | Always Performed |
|---------|-------------------------------|--------------------------------|------------------|
| B190R | BCR/ABL1, p190, Quant, Reflex | Yes, (order BA190) (Bill Only) | No |
| B210R | BCR/ABL1, p210, Quant, Reflex | Yes, (order BCRAW) (Bill Only) | No |

Testing Algorithm

When a positive common p210 or p190 *BCR::ABL1* result is identified by the qualitative assay, a reflex test will then be performed at an additional charge to determine the quantitative transcript level of *BCR::ABL1* messenger RNA. A positive common p210 or p190 result will specifically trigger either quantitative p210 or p190 testing to provide a normalized percentage of transcript level. For the p210 target, the value is additionally defined using the international scale convention. The results are released in an integrated report and provide a baseline quantitative transcript to monitor treatment response. If the initial qualitative testing is negative, or an alternate rare form of *BCR::ABL1* is detected, then no reflex testing will be pursued, and the initial results will be reported.

For more information see:

- [Myeloproliferative Neoplasm: A Diagnostic Approach to Bone Marrow Evaluation](#)
- [Myeloproliferative Neoplasm: A Diagnostic Approach to Peripheral Blood Evaluation](#)

Special Instructions

- [Myeloproliferative Neoplasm: A Diagnostic Approach to Peripheral Blood Evaluation](#)
- [Myeloproliferative Neoplasm: A Diagnostic Approach to Bone Marrow Evaluation](#)
- [Hematopathology Patient Information](#)
- [BCR/ABL1 Ordering Guide for Blood and Bone Marrow](#)

Highlights

Following a positive *BCR::ABL1* diagnostic reverse transcription polymerase chain reaction result, a reflex test will be performed to provide a quantitative measurement of *BCR/ABL1* mRNA transcript (either p190 or p210 types). Current National Comprehensive Cancer Network guidelines for chronic myeloid leukemia, for example, indicate that the quantitative p210 messenger RNA transcript level be obtained at diagnosis. The reflex test establishes the initial patient diagnostic baseline level to assess response to therapy in follow-up samples.

Method Name

Test Definition: BCRFX

BCR/ABL1 Qualitative Diagnostic Assay with
Reflex to BCR/ABL1 p190 Quantitative Assay or
BCR/ABL1 p210 Quantitative Assay, Varies

Reverse Transcription Polymerase Chain Reaction (RT-PCR) Multiplex PCR

NY State Available

Yes

Specimen

Specimen Type

Varies

Ordering Guidance

Additional testing options are available. For ordering guidance see [BCR/ABL1 Ordering Guide for Blood and Bone Marrow](#).

Shipping Instructions

- 1. Specimen must arrive within 72 hours of collection.
- 2. Collect and package specimen as close to shipping time as possible.

Necessary Information

Pertinent clinical history including if the patient has a diagnosis of chronic myeloid/myelogenous leukemia or other *BCR/ABL1* positive neoplasm is required.

Specimen Required

Submit only 1 of the following specimens:

Specimen Type: Whole blood

Container/Tube:

Preferred: Lavender top (EDTA)

Acceptable: Yellow top (ACD)

Specimen Volume: 10 mL

Collection Instructions:

- 1. Invert several times to mix blood.
- 2. Send whole blood specimen in original tube. **Do not aliquot.**
- 3. Label specimen as blood.

Specimen Type: Bone marrow

Container/Tube:

Preferred: Lavender top (EDTA)

Acceptable: Yellow top (ACD)

Specimen Volume: 4 mL

Collection Instructions:

Test Definition: BCRFX

BCR/ABL1 Qualitative Diagnostic Assay with
Reflex to BCR/ABL1 p190 Quantitative Assay or
BCR/ABL1 p210 Quantitative Assay, Varies

1. Invert several times to mix bone marrow.
2. Send bone marrow specimen in original tube. **Do not aliquot.**
3. Label specimen as bone marrow.

Forms

1. [Hematopathology Patient Information](#) (T676)
2. If not ordering electronically, complete, print, and send a [Hematopathology/Cytogenetics Test Request](#) (T726) with the specimen.

Specimen Minimum Volume

Blood: 8 mL; Bone marrow: 2 mL

Reject Due To

| | |
|--------------------------------|--------|
| Gross hemolysis | Reject |
| Moderately to severely clotted | Reject |

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------------|----------|-------------------------|
| Varies | Ambient | 72 hours | PURPLE OR PINK TOP/EDTA |
| | Refrigerated (preferred) | 72 hours | PURPLE OR PINK TOP/EDTA |

Clinical & Interpretive

Clinical Information

The t(9;22)/*BCR::ABL1* abnormality is associated with chronic myelogenous leukemia (CML) and "Philadelphia-positive" acute lymphoblastic leukemia of B-cell lineage (Ph+ ALL). Very rarely, this abnormality has also been identified in cases of acute myeloid leukemia and T-cell lymphoblastic leukemia/lymphoma. The fusion gene on the derivative chromosome 22q11 produces a chimeric *BCR::ABL1* messenger RNA (mRNA) transcript and corresponding translated oncoprotein. Despite substantial breakpoint heterogeneity at the DNA level, a consistent set of *BCR::ABL1* mRNA transcripts are produced that can be readily and sensitively detected by reverse transcription polymerase chain reaction (RT-PCR) technique. In CML, breakpoints in *BCR* result in either exons 13 or 14 (e13, e14) joined to exon 2 of *ABL1* (a2). The corresponding e13-a2 or e14-a2 *BCR::ABL1* mRNAs produce a 210 kDa protein (p210). Rare cases of CML are characterized by an e19-a2 type mRNA with a corresponding p230 protein. In Ph+ ALL, the majority of cases harbor an e1-a2 *BCR::ABL1* mRNA transcript, producing a p190 protein. However, chimeric mRNA type is not invariably associated with disease type, as noted by the presence of p210-positive Ph ALL and very rare cases of p190-positive CML. Therefore, positive results from a screening (diagnostic) assay for *BCR::ABL1* mRNA need to be correlated with clinical and pathologic findings.

In addition to the main transcript variants described above, rare occurrences of both CML and Ph+ ALL can have alternative break-fusion events resulting in unusual *BCR::ABL1* transcript types. Examples include e6-a2 and *BCR* exon fusions to *ABL1* exon a3 (eg, e13-a3, e14-a3, or e1-a3). In addition to detecting common *BCR::ABL1* mRNA transcripts, this assay also can identify these rarer *BCR::ABL1* transcript variants and is therefore a comprehensive screen for both usual and uncommon *BCR::ABL1* gene fusions in hematopoietic malignancies. Given the nature of genetic events in tumors however, this assay will not identify extremely rare and unexpected *BCR::ABL1* events involving other exons (eg, case report level) and is therefore not absolutely specific but is predicted to detect greater than 99.5% of *BCR::ABL1* events. Therefore, it is recommended that for diagnosis, RT-PCR plus a second method (eg, *BCR::ABL1* fluorescence in situ hybridization or cytogenetics) should be used. However, this RT-PCR assay is invaluable at diagnosis for identifying the precise *BCR::ABL1* mRNA type (eg, for future quantitative assay disease monitoring), which complementary methods cannot.

This assay is intended as a qualitative method, providing information on the presence (and specific mRNA type) or absence of the *BCR::ABL1* mRNA. Results from this test can be used to determine the correct subsequent assay for monitoring of transcript levels following therapy (eg, BCRA1 / *BCR/ABL1*, p210, mRNA Detection, Reverse Transcription-PCR (RT-PCR), Quantitative, Monitoring Chronic Myeloid Leukemia (CML), Varies; BA190 / *BCR/ABL1*, p190, mRNA Detection, Reverse Transcription-PCR (RT-PCR), Quantitative, Monitoring Assay, Varies). Because the assay is analytically sensitive, it compensates for situations such as partially degraded RNA quality or low cell number, but it is not intended for quantitative or monitoring purposes.

Reference Values

An interpretive report will be provided.

Interpretation

An interpretive report will be provided.

When positive, the test identifies which specific messenger RNA fusion variant is present to guide selection of an appropriate monitoring assay. If common p210 or p190 fusion variant detected, quantitative reflex will be performed.

-Common fusion variants detected: e13-a2 or e14-a2 (p210), e1-a2 (p190), and e6-a2 (p205*)

-Rare fusion variants detected: e13-a3 (p210), e14-a3 (p210), e1-a3 (p190), e19-a2 (p230)

-Potential rare fusions detected: e12-a3, e19-a3

*This is formerly observed as the e6-a2 (p185) fusion form

Cautions

No significant cautionary statements

Clinical Reference

1. Burmeister T, Reinhardt R. A multiplex PCR for improved detection of typical and atypical BCR-ABL fusion transcripts. *Leuk Res*. 2008;32(4):579-585
2. Melo JV. The diversity of BCR-ABL fusion proteins and their relationship to leukemia phenotype. *Blood*. 1996;88(7):2375-2384
3. Melo JV. BCR-ABL gene variants. *Baillieres Clin Haematol*. 1997;10(2):203-222
4. Baccarini M, Deininger MW, Rosti G, et al. European LeukemiaNet recommendations for the management of chronic

myeloid leukemia: 2013. Blood. 2013;122(6):872-884

5. Cross NC, White HE, Muller MC, et al. Standardized definitions of molecular response in chronic myeloid leukemia. Leukemia. 2012;26(10):2172-2175

6. Deininger MW, Shah NP, Altman JK, et al. Chronic Myeloid Leukemia, Version 2.2021, NCCN Clinical Practice Guidelines in Oncology. J Natl Compr Canc Netw. 2020;18(10):1385-1415. doi:10.6004/jnccn.2020.0047

Performance

Method Description

Total RNA is extracted from the patient's blood or bone marrow at the time of diagnosis and messenger RNA (mRNA) is reverse transcribed into complementary DNA (cDNA). The cDNA is then subjected to polymerase chain reaction (PCR) using 4 separate multiplex reactions. A qualitative result, which will include the relative ratio of target translocation mRNA to control *GUSB* gene mRNA, will be given by LightCycler 96 software. As this method employs a quantitative PCR platform, the results can be used to evaluate the relative expression levels of the translocation mRNA relative to control mRNA, thus providing an improved measure of RNA quality in the assay. Reporting of results will be qualitative; either *BCR::ABL1* mRNA positive/detected (with transcript type) or negative/not detected.(Unpublished Mayo method)

PDF Report

No

Day(s) Performed

Monday through Saturday

Report Available

7 to 10 days

Specimen Retention Time

Blood, Bone marrow: 2 weeks; Extracted RNA: 3 months

Performing Laboratory Location

Rochester

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

Test Definition: BCRFX

BCR/ABL1 Qualitative Diagnostic Assay with
Reflex to BCR/ABL1 p190 Quantitative Assay or
BCR/ABL1 p210 Quantitative Assay, Varies

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

CPT Code Information

- 81206
- 81207
- 81208
- 81479 (if appropriate for government payers)

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|-----------------------------|--------------------|
| BCRFX | BCR/ABL1 Reflex, Qual/Quant | In Process |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------------|---------------------|
| MP039 | Specimen Type | 31208-2 |
| 48389 | BCR/ABL1 Reflex Result | No LOINC Needed |
| 48388 | Interpretation | 69047-9 |