

Overview

Useful For

Diagnosing vitamin A deficiency and toxicity

Monitoring vitamin A therapy

Method Name

Liquid Chromatography-Tandem Mass Spectrometry (LC-MS/MS)

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Patient Preparation:

**Fasting:** 12 hours, required; infants should have specimen collected before next feeding.

**Supplies:** Sarstedt Aliquot Tube, 5 mL (T914)

Collection Container/Tube:

**Preferred:** Red top

**Acceptable:** Serum gel

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 0.5 mL

**Collection Instructions:** Within 2 hours of collection, centrifuge and aliquot serum into a plastic vial.

Note: Serum must be removed from the red blood cells within 2 hours of collection.

Forms

If not ordering electronically, complete, print, and send a [General Request](#) (T239) with the specimen.

Specimen Minimum Volume

0.25 mL

Reject Due To

Gross hemolysis	OK
Gross lipemia	OK
Gross icterus	OK

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	28 days	

Clinical & Interpretive

Clinical Information

The level of vitamin A in the plasma or serum is a reflection of the quantities of vitamin A and carotene (provitamin A) ingested and absorbed by the intestine (carotene is converted to vitamin A by intestinal absorptive cells and hepatocytes).

Vitamin A plays an essential role in the function of the retina (adaptation to dim light), is necessary for growth and differentiation of epithelial tissue, and is required for growth of bone, reproduction, and embryonic development. Together with certain carotenoids, vitamin A enhances immune function, reducing the consequences of some infectious diseases.

Degenerative changes in eyes and skin are commonly observed in vitamin A deficiency. Poor adaptation of vision to darkness (night blindness) is an early symptom that may be followed by degenerative changes in the retina. In developing countries, vitamin A deficiency is the principal preventable cause of blindness. Severe or prolonged deficiency leads to dry eye (xerophthalmia) that can result in corneal ulcers, scarring, and blindness. Another important consequence of inadequate intake is acquired immunodeficiency disease, where an increased incidence of death is associated with deficient vitamin A levels. Increased susceptibility is associated with vitamin A deficiency. In patients with HIV, vitamin A deficiency is associated with increased disease progression and mortality.

Vitamin A in excess can be toxic. In particular, chronic vitamin A intoxication is a concern in normal adults who ingest more than 15 mg per day and children who ingest more than 6 mg per day of vitamin A over a period of several months. Manifestations are various and include dry skin, cheilosis, glossitis, vomiting, alopecia, bone demineralization and pain, hypercalcemia, lymph node enlargement, hyperlipidemia, amenorrhea, and features of pseudotumor cerebri with increased intracranial pressure and papilledema. Liver fibrosis with portal hypertension may also result. Congenital malformations, like spontaneous abortions, craniofacial abnormalities, and valvular heart disease have been described in pregnant women taking vitamin A in excess. Consequently, in pregnancy, the daily dose of vitamin A should not exceed 3 mg.

Reference Values

- 0-6 years: 11.3-64.7 mcg/dL
- 7-12 years: 12.8-81.2 mcg/dL
- 13-17 years: 14.4-97.7 mcg/dL
- > or =18 years: 32.5-78.0 mcg/dL

Interpretation

The World Health Organization recommends supplementation when vitamin A levels fall below 20.0 mcg/dL.

Severe deficiency is indicated at levels less than 10.0 mcg/dL.

Vitamin A values above 120.0 mcg/dL suggest hypervitaminosis A and associated toxicity.

### Cautions

Acute ethanol ingestion may result in increased serum vitamin A levels.

Testing of nonfasting specimens or the use of vitamin supplementation can result in elevated serum vitamin concentrations. Reference values were established using specimens from individuals who were fasting.

### Clinical Reference

1. Ball GFM. Vitamins: Their role in the human body. Blackwell Publishing; 2004:133-187
2. Ross AC. Vitamin A and carotenoids. In: Shils ME, Shike M, Ross AC, et al, eds. Modern Nutrition in Health and Disease. 10th ed. Lippincott Williams and Wilkins; 2006:351-375
3. Roberts NB, Taylor A, Sodi R. Vitamins and Trace Elements. In: Rifai N, Horvath AR, Wittwer CT, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 6th ed. Elsevier; 2018:chap 37
4. Sodi R. Vitamins and trace elements. In: Rifai N, Chiu RWK, Young I, Burnham C-AD, Wittwer CT, eds. Tietz Textbook of Laboratory Medicine. 7th ed. Elsevier; 2023:417-417.e104

## Performance

### Method Description

Deuterated vitamin A (d6-all-trans retinol) is added to serum as an internal standard. Vitamin A (all-trans retinol) and the deuterated internal standard are extracted from the specimens using on-line turbulent flow high-performance liquid chromatography and analyzed by liquid chromatography-tandem mass spectrometry using multiple reaction monitoring in positive mode. (Unpublished Mayo method)

### PDF Report

No

### Day(s) Performed

Sunday through Friday

### Report Available

2 to 5 days

### Specimen Retention Time

14 days

### Performing Laboratory Location

Mayo Clinic Laboratories - Rochester Superior Drive

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

CPT Code Information

84590

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
VITA	Vitamin A, S	2923-1

Result ID	Test Result Name	Result LOINC® Value
7597	Vitamin A	2923-1