

## Overview

### Method Name

Quantitative Enzymatic

### NY State Available

Yes

## Specimen

### Specimen Type

Serum

### Specimen Required

**Specimen Type:** Serum

**Collection Container/Tube:** SST or Red

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

#### Collection Instructions:

1. Draw blood in a serum gel tube(s), plain red-top tube(s) is acceptable.
2. Centrifuge and send 1 mL of serum refrigerate in a plastic vial.

### Specimen Minimum Volume

0.5 mL

### Reject Due To

Hemolysis:	Mild; reject, Gross; reject
Thawing:	Warm reject; Cold OK
Lipemia:	NA
Icterus:	NA
Other:	NA

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	14 days	
	Ambient	4 hours	

**Clinical & Interpretive****Reference Values**

0 - 15 U/L

**Performance****PDF Report**

No

**Day(s) Performed**

Sunday through Saturday

**Report Available**

1 to 6 days

**Performing Laboratory Location**

ARUP Laboratories

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by ARUP Laboratories. The U.S. Food and Drug Administration has not approved or cleared this test; however, FDA clearance or approval is not currently required for clinical use. The results are not intended to be used as the sole means for clinical diagnosis or patient management decisions.

**CPT Code Information**

83915

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
F5NUL	5' Nucleotidase	1690-7

Result ID	Test Result Name	Result LOINC® Value
Z2607	5' Nucleotidase	1690-7