

**Overview****Method Name**

Quantitative Enzyme-Linked Immunosorbent Assay

**NY State Available**

Yes

**Specimen****Specimen Type**

WB Heparin

**Shipping Instructions**

CRITICAL FROZEN. Separate specimens must be submitted when multiple tests are ordered.

**Specimen Required****Specimen Type:** Whole Blood**Collection Container/Tube:****Preferred:** Green top (sodium heparin)**Acceptable:** Green top (lithium heparin)**Submission Container/Tube:** Plastic vial**Specimen Volume:** 1 mL**Collection Instructions:**

1. Draw blood in a green-top tube (sodium or lithium heparin).
2. Aliquot 1 mL well mixed whole blood into a plastic vial.
3. Freeze specimen.

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

Ambient or refrigerated samples	Reject
---------------------------------	--------

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
WB Heparin	Frozen	180 days	

**Clinical & Interpretive****Reference Values**

180 - 1800 nmol/L

**Performance****PDF Report**

No

**Day(s) Performed**

Monday, Thursday

**Report Available**

3 to 11 days

**Performing Laboratory Location**

ARUP Laboratories

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

**CPT Code Information**

83088

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FHSTW	Histamine, Whole Blood	46436-2

Result ID	Test Result Name	Result LOINC® Value
-----------	------------------	---------------------

---

Z2753	Histamine, Whole Blood	46436-2
-------	------------------------	---------