

Test Definition: FANSE

Annatto Seed (Bixa orellana) IgE

Overview

Method Name

Enzyme Immunoassay (EIA)

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Container/Tube: Red or SST Preferred: Red top tube Acceptable: Serum gel tube Specimen Volume: 1 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5

mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.5 mL

Reject Due To

Hemolysis	NA NA
Lipemia	Mild OK; Gross reject
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

Clinical & Interpretive



Test Definition: FANSE

Annatto Seed (Bixa orellana) IgE

Reference Values

<0.35 kU/L

Interpretation

Class	IgE (kU/L)	<u>Comment</u>
0	<0.35	Below Detection
1	0.35 -0.69	Low Positive
2	0.70 - 3.49	Moderate Positive
3	3.50 – 17.49	Positive
4	17.50 – 49.99	Strong Positive
5	50.00 – 99.99	Very Strong Positive
6	>99.99	Very Strong Positive

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

2 to 7 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.



Test Definition: FANSE

Annatto Seed (Bixa orellana) IgE

CPT Code Information

86003

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FANSE	Annatto Seed IgE	50638-6

Result ID	Test Result Name	Result LOINC® Value	
Z3074	Annatto Seed IgE	50638-6	
Z3075	Class	103035-2	