

# **Test Definition: FPOPW**

Poplar White (Populus alba) IgE

Overview	
	١

## **Method Name**

Enzyme Immunoassay (EIA)

#### **NY State Available**

Yes

# **Specimen**

# **Specimen Type**

Serum

## **Specimen Required**

Draw blood in a plain, red-top tube, serum gel tube are acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.

# **Specimen Minimum Volume**

0.5 mL

## **Reject Due To**

Thawing:	Warm OK; Cold OK
Lipemia:	Reject

## **Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

# **Clinical & Interpretive**

## **Reference Values**

<0.35 kU/L

## Interpretation



# **Test Definition: FPOPW**

Poplar White (Populus alba) IgE

<u>Class</u>	IgE (kU/L)	Comment
0	<0.35	Below Detection
1	0.35-0.69	Low Positive
2	0.70-3.49	Moderate Positive
3	3.50-17.49	Positive
4	17.50-49.99	Strong Positive
5	>49.99	Very Strong Positive

#### **Performance**

#### **PDF Report**

No

### Day(s) Performed

Monday through Friday

#### **Report Available**

4 to 7 days

#### **Performing Laboratory Location**

**Eurofins Viracor** 

#### **Fees & Codes**

#### **Fees**

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

#### **Test Classification**

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

#### **CPT Code Information**

86003

# **LOINC®** Information

Test ID	Test Order Name	Order LOINC® Value
FPOPW	Poplar White (Populus alba) IgE	21473-4
Result ID	Test Result Name	Result LOINC® Value



# **Test Definition: FPOPW**

Poplar White (Populus alba) IgE

Z3114	Poplar White (Populus alba) IgE	21473-4
Z3115	Class	102802-6