

**Overview****Method Name**

ImmunoCAP FEIA

**NY State Available**

Yes

**Specimen****Specimen Type**

Serum

**Specimen Required****Specimen Type:** Serum**Container/Tube:** Red Top**Specimen Volume:** 0.5 mL**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

Thawing:	Warm OK; Cold OK
Lipemia:	Reject

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Frozen	365 days	
	Ambient	28 days	

**Clinical & Interpretive****Reference Values**

---

<0.35 kU/L

**Interpretation**

<a href="#">Class</a>	IgE (kU/L)	Comment
0	<0.10	Negative
0/1	0.10-0.34	Equivocal
1	0.35-0.69	Low Positive
2	0.70-3.49	Moderate Positive
3	3.50-17.49	High Positive
4	17.50-49.99	Very High Positive
5	50.00-99.99	Very High Positive
6	>99.99	Very High Positive

**Performance****PDF Report**

No

**Day(s) Performed**

Monday through Friday

**Report Available**

3 to 6 days

**Performing Laboratory Location**

Eurofins Viracor

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

[This test was developed and its performance characteristics determined by Eurofins Viracor . It has not been cleared or approved by the U.S. Food and Drug Administration.](#)

**CPT Code Information**

## Test Definition: FWALP

Wall Eyed Pike (Sander vitreus)(Stizostedium  
vitreum) IgE

---

86003

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FWALP	Wall-Eyed Pike IgE	7761-0

Result ID	Test Result Name	Result LOINC® Value
Z3120	Wall-Eyed Pike IgE	7761-0
Z3121	Class	102796-0