

# **Test Definition: FBWME**

Bayberry/Wax Myrtle (Myrica spp) IgE

## **Overview**

#### **Method Name**

ImmunoCAP (FEIA)

#### **NY State Available**

Yes

# **Specimen**

## **Specimen Type**

Serum

# **Specimen Required**

Specimen Type: Serum
Container/Tube: Red Top
Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain re-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL

serum refrigerated in a plastic vial.

#### **Specimen Minimum Volume**

0.5 mL

#### **Reject Due To**

Lipemia	Reject

## **Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

# **Clinical & Interpretive**

#### **Reference Values**

<0.35 kU/L



# **Test Definition: FBWME**

Bayberry/Wax Myrtle (Myrica spp) IgE

#### Interpretation

<u>Class</u>	IgE (kU/L)	<u>Comment</u>
0	<0.10	Negative
0/1	0.10 - 0.34	Equivocal
1	0.35 - 0.69	Low Positive
2	0.70 - 3.49	Moderate Positive
3	3.50 - 17.49	High Positive
4	17.50 - 49.99	Very High Positive
5	50.00 - 99.99	Very High Positive
6	>99.99	Very High Positive

#### **Performance**

#### **PDF Report**

No

#### Day(s) Performed

Monday through Friday

#### **Report Available**

3 to 6 days

#### **Performing Laboratory Location**

Eurofins Viracor

#### **Fees & Codes**

#### **Fees**

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

#### **Test Classification**

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

#### **CPT Code Information**

86003

#### **LOINC®** Information



# **Test Definition: FBWME**

Bayberry/Wax Myrtle (Myrica spp) IgE

Test ID	Test Order Name	Order LOINC® Value
FBWME	Bayberry/Wax Myrtle IgE	7123-3

Result ID	Test Result Name	Result LOINC® Value
Z3130	Bayberry/Wax Myrtle IgE	7123-3
Z3131	Class	102164-1