

Overview

Method Name

ImmunoCAP (FEIA)

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

**Specimen Type:** Serum

**Container/Tube:** Red Top

**Specimen Volume:** 0.5 mL

**Collection Instructions:** Draw blood in a plain re-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.5 mL

Reject Due To

Lipemia	Reject
---------	--------

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

Clinical & Interpretive

Reference Values

<0.35 kU/L

Interpretation

Class	IgE (kU/L)	Comment
0	<0.10	Negative
0/1	0.10 - 0.34	Equivocal
1	0.35 - 0.69	Low Positive
2	0.70 - 3.49	Moderate Positive
3	3.50 - 17.49	High Positive
4	17.50 - 49.99	Very High Positive
5	50.00 - 99.99	Very High Positive
6	>99.99	Very High Positive

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

3 to 6 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

86003

LOINC® Information

# Test Definition: FBWME

Bayberry/Wax Myrtle (Myrica spp) IgE

Test ID	Test Order Name	Order LOINC® Value
FBWME	Bayberry/Wax Myrtle IgE	7123-3

Result ID	Test Result Name	Result LOINC® Value
Z3130	Bayberry/Wax Myrtle IgE	7123-3
Z3131	Class	102164-1