

Overview

Method Name

Qualitative Hemagglutination

NY State Available

Yes

Specimen

Specimen Type

Plasma Na Cit

Specimen Required

Collection Container/Tube: Light-blue top (3.2% sodium citrate)

Submission Container/Tube: Plastic vial

Specimen Volume: 1.5 mL Platelet-poor plasma

Collection Instructions:

1. Centrifuge at 1700 x *g* for 15 minutes.
2. Aliquot 1.5 platelet-poor plasma into plastic vial making sure not to disturb the buffy coat.
3. Send frozen.

Note: Critical frozen. Separate specimens must be submitted when multiple tests are ordered.

Specimen Minimum Volume

Platelet-poor plasma: 1 mL

Reject Due To

| | |
|--------------------|--------|
| Hemolysis | Reject |
| Clotted | Reject |
| Glass vials | Reject |
| Serum, EDTA plasma | Reject |

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|-------------|---------|-------------------|
| Plasma Na Cit | Frozen | 14 days | |

Clinical & Interpretive

Clinical Information

Refer to <https://ltd.aruplab.com/>

Reference Values

Negative

Cautions

Activation of specimen due to poor venipuncture technique may cause false-positive results.

Performance**PDF Report**

No

Day(s) Performed

Monday through Sunday

Report Available

1 to 6 days

Performing Laboratory Location

ARUP Laboratories

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

85366

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|------------------------|--------------------|
| FSFM | Soluble Fibrin Monomer | 40702-3 |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------------|---------------------|
| FSFM | Soluble Fibrin Monomer | 40702-3 |