

## Overview

### Useful For

Identifying solid tumors that may respond to targeted therapies by simultaneously assessing for fusions involving targeted regions of the *NTRK1*, *NTRK2*, and *NTRK3* genes resulting in fusion transcripts.

This test is **not useful for** hematologic malignancies.

### Genetics Test Information

This test uses next-generation sequencing to identify rearrangements (fusions) involving targeted regions of the *NTRK1*, *NTRK2*, and *NTRK3* genes.

The assay was designed to detect the most common *NTRK* gene fusions and novel gene fusion partners in the targeted gene regions. For more information see [NTRK RNA Targeted Gene Fusions](#).

Of note, this test is performed to evaluate rearrangements (fusions) within solid tumor samples. This test is not intended for use for hematological malignancies.

### Additional Tests

Test Id	Reporting Name	Available Separately	Always Performed
SLIRV	Slide Review in MG	No, (Bill Only)	Yes

### Testing Algorithm

When this test is ordered, slide review will always be performed at an additional charge.

### Special Instructions

- [NTRK RNA Targeted Gene Fusions](#)
- [Tissue Requirements for Solid Tumor Next-Generation Sequencing](#)

### Highlights

This test evaluates formalin-fixed, paraffin-embedded tumor slides from patients with advanced solid tumors for rearrangements (fusions) involving targeted regions of the *NTRK1*, *NTRK2*, and *NTRK3* genes to identify candidates for targeted therapy. Current data suggests that solid tumors with *NTRK* rearrangements may be sensitive to multikinase inhibitors.

### Method Name

Polymerase Chain Reaction (PCR)-based Next-Generation Sequencing (NGS)

### NY State Available

Yes

## Specimen

### Specimen Type

Varies

### Ordering Guidance

Multiple oncology (cancer) gene panels are available. For more information see [Hematology, Oncology, and Hereditary Test Selection Guide](#).

### Necessary Information

Pathology report (final or preliminary) at minimum containing the following information must accompany specimen in order for testing to be performed:

1. Patient name
2. Block number-must be on all blocks, slides and paperwork (can be handwritten on the paperwork)
3. Tissue collection date
4. Source of the tissue

### Specimen Required

**This assay requires at least 10% tumor nuclei.**

-Preferred amount of tumor area with sufficient percent tumor nuclei: tissue 144 mm<sup>2</sup>)

-Minimum amount of tumor area: tissue 36 mm<sup>2</sup>).

-These amounts are cumulative over up to 10 unstained slides and must have adequate percent tumor nuclei.

-Tissue fixation: 10% neutral buffered formalin, not decalcified

-For specimen preparation guidance, see [Tissue Requirement for Solid Tumor Next-Generation Sequencing](#). In this document, the sizes are given as 4 mm x 4 mm x 10 slides as preferred: approximate/equivalent to 144 mm<sup>2</sup>) and the minimum as 3 mm x 1 mm x 10 slides: approximate/equivalent to 36 mm<sup>2</sup>).

#### Preferred:

**Specimen Type:** Tissue block

**Collection Instructions:** Submit a formalin-fixed, paraffin-embedded tissue block with acceptable amount of tumor tissue.

#### Acceptable:

**Specimen Type:** Tissue slide

**Slides:** 1 stained and 10 unstained

**Collection Instructions:** Submit 1 slide stained with hematoxylin and eosin and 10 unstained, nonbaked slides with 5-micron thick sections of the tumor tissue.

**Note:** The total amount of required tumor nuclei can be obtained by scraping up to 10 slides from the same block.

### Forms

If not ordering electronically, complete, print, and send a [Oncology Test Request](#) (T729) with the specimen.

### Specimen Minimum Volume

See Specimen Required

**Reject Due To**

Specimens that have been decalcified (all methods) Specimens that have not been formalin-fixed, paraffin-embedded Bone marrow in EDTA	Reject
---	--------

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)		
	Refrigerated		
	Frozen		

**Clinical & Interpretive**

**Clinical Information**

Targeted cancer therapies are defined as antibody or small molecule drugs that block the growth and spread of cancer by interfering with specific cell molecules involved in tumor growth and progression. Multiple targeted therapies have been approved by the US Food and Drug Administration for treatment of specific cancers. Molecular genetic profiling is often needed to identify targets amenable to targeted therapies and to minimize treatment costs and therapy-associated risks.

Fusions involving the *NTRK1*, *NTRK2*, and *NTRK3* genes (ie, NTRK gene fusions) are oncogenic drivers of multiple types of pediatric and adult solid tumors. In solid tumors, the presence of an *NTRK* gene fusion is a biomarker for response to tropomyosin receptor kinase (TRK) inhibitor therapies.

This test assesses for fusions involving targeted regions of *NTRK1*, *NTRK2*, and *NTRK3* genes. The results of this test can be useful in guiding treatment of individuals with advanced solid tumors.

See Method Description for details regarding the targeted gene regions evaluated by this test.

**Reference Values**

An interpretive report will be provided.

**Interpretation**

---

The interpretation of molecular biomarker analysis includes an overview of the results and the associated diagnostic, prognostic, and therapeutic implications.

**Cautions**

This test is not designed to differentiate between somatic and germline alterations. Additional testing may be necessary to clarify the significance of results if there is a potential hereditary risk.

A negative (wildtype) result does not rule out the presence of a rearrangement (fusion) that may be present but below the limits of detection of this assay. The analytical sensitivity of this assay is a minimum coverage of 10 targeted fusion reads with 5 unique fusion molecules in a sample with 10% or greater tumor content.

Only gene rearrangements (fusions) involving targeted regions in *NTRK1*, *NTRK2*, and *NTRK3* genes will be detected. This test does not detect point mutations, insertion/deletion mutations, large single or multiexon deletions or duplications, or genomic copy number variants in any of the genes tested.

Rare alterations (ie, polymorphisms) may be present that could lead to false-negative or false-positive results.

Test results should be interpreted in the context of clinical findings, tumor sampling, and other laboratory data. If results obtained do not match other clinical or laboratory findings, contact the laboratory for an updated interpretation. Misinterpretation of results may occur if the information provided is inaccurate or incomplete.

Reliable results are dependent on adequate specimen collection and processing. This test has been validated on formalin-fixed, paraffin-embedded tissues; other types of fixatives are discouraged. Improper treatment of tissues, such as decalcification, may cause polymerase chain reaction failure.

**Clinical Reference**

1. Drilon A, Laetsch TW, Kummar S. Efficacy of Larotrectinib in TRK fusion-positive cancers in adults and children. *N Engl J Med*. 2018;378(8):731-739. doi:10.1056/NEJMoa1714448
2. Cocco E, Scaltriti M, Drilon A. NTRK fusion-positive cancers and TRK inhibitor therapy. *Nat Rev Clin Oncol*. 2018;15(12):731-747. doi:10.1038/s41571-018-0113-0

**Performance****Method Description**

Next-generation sequencing (NGS) is performed to test for the presence of fusions involving targeted regions in *NTRK1*, *NTRK2*, and *NTRK3* genes. Fusion gene nomenclature is based on a custom reference sequence using genome build GRCh37 (hg19). See [NTRK RNA Targeted Gene Fusions](#) for details regarding the targeted gene regions identified by this test. (Unpublished Mayo method)

**PDF Report**

No

**Day(s) Performed**

Varies

**Report Available**

14 to 21 days

**Specimen Retention Time**

FPPE tissue block: Unused portions of FPPE blocks will be returned; Unused, unstained slides: 5 years; Stained slides: Indefinitely

**Performing Laboratory Location**

Rochester

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

**CPT Code Information**

81194

88381

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
NTRK	NTRK Gene Fusion Panel	93813-4

Result ID	Test Result Name	Result LOINC® Value
606761	Result Summary	50397-9
606762	Result	81286-7
606763	Interpretation	69047-9
606764	Additional Information	48767-8
606765	Specimen	31208-2
606766	Source	31208-2
606767	Tissue ID	80398-1
606768	Released By	18771-6