

---

## Overview

### Useful For

Genetic test for individuals at high risk for glucose-6-phosphate dehydrogenase (G6PD) deficiency

Aiding in the diagnosis of G6PD deficiency

Determining G6PD deficiency status in individuals with inconclusive or unexpected phenotyping results

Differentiation of heterozygotes with skewed X-inactivation from homozygotes and compound heterozygotes

Definitive diagnosis of carrier status

Evaluation of neonates with unexplained jaundice

Identifying individuals at risk of drug-induced acute hemolytic anemia related to G6PD deficiency

### Genetics Test Information

**This test is for molecular sequencing of the *G6PD* gene and does not assess glucose-6-phosphate dehydrogenase (G6PD) enzyme activity.** Enzymatic testing may be suggested as follow-up to this assay. For G6PD enzyme testing order G6PD1 / Glucose 6-Phosphate Dehydrogenase Enzyme Activity, Blood.

G6PD deficiency is a common X-linked condition, estimated to affect up to 500 million people worldwide. Both male and female patients may be impacted due to how common G6PD deficiency is in the population.

Acute hemolytic anemia (AHA) can be triggered in individuals with G6PD deficiency by fava beans, several types of medications (including rasburicase, dapsone-containing combinations of antimalarial drugs, and methylene blue), and infection. Less commonly, chronic congenital nonspherocytic hemolytic anemia (CNSHA) may occur in severe forms of G6PD deficiency.

US Food and Drug Administration labeling and Clinical Pharmacogenetics Implementation Consortium (CPIC) guidelines recommend that G6PD testing be undertaken in high-risk populations before prescribing drugs known to cause AHA. Knowing a patient's genotype is generally sufficient to avoid contraindicated drugs, but follow-up with the phenotyping (enzyme) assay may be necessary to clarify results in some cases.

This test involves full gene sequencing of all exons and exon/intron boundaries of the *G6PD* gene. A comprehensive interpretation will be provided including congenital and pharmacogenomic implications of results. Testing should be considered before prescribing medication associated with hemolysis in individuals with G6PD deficiency.

### Testing Algorithm

The following are available:

[-Glucose-6-Phosphate Dehydrogenase \(G6PD\) Deficiency Diagnostic Algorithm](#)

---

[-Glucose-6-Phosphate Dehydrogenase \(G6PD\) Genotyping Interpretive Algorithm](#)

**Special Instructions**

- [Informed Consent for Genetic Testing](#)
- [Pharmacogenomic Association Tables](#)
- [Glucose-6-Phosphate Dehydrogenase \(G6PD\) Genotyping Interpretive Algorithm](#)
- [Glucose-6-Phosphate Dehydrogenase \(G6PD\) Deficiency Diagnostic Algorithm](#)

**Method Name**

Polymerase Chain Reaction (PCR) followed by DNA Sequence Analysis

**NY State Available**

Yes

**Specimen****Specimen Type**

Varies

**Ordering Guidance**

For initial or time-sensitive screening for glucose-6-phosphate dehydrogenase deficiency, order G6PD1 / Glucose 6-Phosphate Dehydrogenase Enzyme Activity, Blood.

**Necessary Information**

Include physician name and phone number with the specimen.

**Specimen Required**

Submit only 1 of the following specimens:

**Specimen Type:** Whole blood

**Container/Tube:**

**Preferred:** Lavender top (EDTA) or yellow top (ACD)

**Specimen Volume:** 3 mL

**Collection Instructions:**

1. Invert several times to mix blood.
2. Send whole blood specimen in original tube. **Do not** aliquot.

**Specimen Stability Information:** Ambient (preferred) 9 days/Refrigerated 30 days

**Specimen Type:** Saliva

**Patient Preparation:** Patient should not eat, drink, smoke, or chew gum 30 minutes prior to collection.

**Supplies:** Saliva Swab Collection Kit (T786)

**Specimen Volume:** 1 Swab

**Collection Instructions:** Collect and send specimen per kit instructions.

**Specimen Stability Information:** Ambient 30 days

### Forms

1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file. The following documents are available:

-[Informed Consent for Genetic Testing](#) (T576)

-[Informed Consent for Genetic Testing \(Spanish\)](#) (T826)

2. If not ordering electronically, complete, print, and send 1 of the following forms with the specimen:

-[Therapeutics Test Request](#) (T831)

-[Benign Hematology Test Request Form](#) (T755)

### Specimen Minimum Volume

Blood: 0.45 mL

Saliva: See Specimen Required

### Reject Due To

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Varies		

### Clinical & Interpretive

#### Clinical Information

Glucose-6-phosphate dehydrogenase (G6PD) deficiency is the most common human enzymopathy, estimated to affect up to 500 million people worldwide. It is most frequently found in populations where *Plasmodium falciparum* malaria is (or was) endemic, but G6PD deficiency may be present in any population.

G6PD deficiency primarily manifests as episodic acute hemolytic anemia (AHA), chronic non-spherocytic hemolytic anemia (CNSHA), and neonatal jaundice. These clinical manifestations can be triggered in individuals with G6PD deficiency by fava beans, several types of medications (including rasburicase, dapsone-containing combinations of antimalarial drugs, and methylene blue), and infection.

G6PD converts glucose-6-phosphate to 6-phosphoglyconolactone in the first step of the pentose phosphate pathway, this reaction also produces nicotinamide adenine dinucleotide phosphate (NADPH) from NADP(+). NADPH, through subsequent enzymatic reactions, protects erythrocytes from damage by detoxifying hydrogen peroxide and other sources of oxidative stress.

G6PD is encoded by the gene *G6PD*, which lies on the X-chromosome. G6PD deficiency is inherited in an X-linked

recessive manner; therefore, male patients are more commonly affected than female patients, but due to the high prevalence of G6PD deficiency, homozygous and compound heterozygous female patients are not uncommon. Over 200 *G6PD* variants have been discovered and are classified based on guidance from the World Health Organization (WHO). In 2022, WHO proposed updated guidance for the classification of *G6PD* variants (Table). This revised guidance is based on the median residual enzyme activity and seeks to resolve problems identified with the WHO *G6PD* classification system that has been in place since 1985 (Table).

Table. Updated and Legacy *G6PD* Variant WHO Classification and Associated G6PD Deficiency Phenotype

2022 WHO class	Median G6PD activity	Hemolysis	Legacy WHO class	Level of residual enzyme activity (% of normal)
A	<20%	Chronic (CNSHA)	I	<10%
B	<45%	Acute, triggered	II	<10%
			III	10%-60%
C	60-150%	No hemolysis	IV	Normal
U	Any	Uncertain clinical significance		

With the exception of those with CNSHA, individuals with G6PD deficiency are typically asymptomatic until they are challenged with an exogenous factor, such as a drug, infection, or fava beans. The exogenous factor can trigger AHA in individuals with G6PD deficiency. The severity of AHA is highly variable, ranging from mild neonatal jaundice to life-threatening complications, such as kernicterus. Therefore, determining the G6PD deficiency status is recommended on the US Food and Drug Administration label of several drugs either proven or suspected to cause AHA in patients with G6PD deficiency. For more information on drugs known to cause AHA in individuals with G6PD deficiency, see [Pharmacogenomic Associations Tables](#).

Preemptive genotyping allows for the identification of patients at risk for an adverse reaction to drugs known to cause AHA in those with G6PD deficiency. In most cases, genotyping provides sufficient information to avoid the use of contraindicated drugs. In some cases, including heterozygous female patients, the phenotyping assay is necessary to determine if such drugs should be avoided. Skewed X-inactivation in heterozygous female patients has been reported to result in G6PD deficiency, so the phenotyping assay is necessary to determine G6PD activity level.

### Reference Values

An interpretive report will be provided.

### Interpretation

All detected alterations will be evaluated according to the latest American College of Medical Genetics and Genomics recommendations and the most recent World Health Organization system for classifying genetic variants of *G6PD*.<sup>(1,2)</sup> Variants will be classified based on known, predicted, or possible effect on gene pathogenicity and reported with interpretive comments detailing their potential or known significance.

### Cautions

Patients who have received a non-leukocyte reduced blood transfusion within the preceding 6 weeks, or who have

---

received an allogeneic hematopoietic stem cell transplant, can have inaccurate genetic test results due to the presence of both donor and recipient DNA.

For patients who have been transfused within the preceding 6 weeks, the glucose-6-phosphate dehydrogenase (G6PD) enzyme assay will also be affected, so it is **not** an appropriate alternative test.

Patients who have received an allogeneic hematopoietic stem cell transplant would be expected to convert G6PD status to that of donor. However, if the patient's transplant was partially successful or if there is a relapse of an underlying hematologic malignancy, a mixture of donor and recipient genotype may be seen on genetic analysis. The enzyme assay can be run after transplantation; order G6PD1 / Glucose 6-Phosphate Dehydrogenase Enzyme Activity, Blood.

Rare variants exist that could lead to false-negative or false-positive results. Other variants in the primer binding regions can affect the testing, and ultimately, the genotype assessment made.

Test results should be interpreted in the context of clinical findings, family history, and other laboratory data. Large deletions or rearrangements are not detected by this assay.

Sometimes a genetic alteration of unknown significance may be identified. In this case, testing of appropriate family members may be useful to determine pathogenicity of the alteration.

This test is not designed to provide specific dosing or drug selection recommendations and is to be used as an aid to clinical decision making only. Drug-label guidance should be used when dosing patients with medications regardless of the predicted phenotype.

Skewed X-inactivation in heterozygous female patients has been reported to result in G6PD deficiency. In these cases, the phenotyping (enzyme) assay is necessary to determine G6PD activity level and assign G6PD deficiency status.

Rarely, incidental or secondary findings may implicate another predisposition or presence of active disease. Incidental findings may include, but are not limited to, results related to the sex chromosomes. These findings will be carefully reviewed to determine whether they will be reported.

### Clinical Reference

1. Richards S, Aziz N, Bale S, et al. Standards and guidelines for the interpretation of sequence variants: a joint consensus recommendation of the American College of Medical Genetics and Genomics and the Association for Molecular Pathology. *Genet Med*. 2015;(17):105-423
2. Global Malaria Programme, Malaria Policy Advisory Group. Meeting report of the technical consultation to review the classification of glucose-6-phosphate dehydrogenase (G6PD). World Health Organization; 2022. Accessed July 3, 2023. Available at [www.who.int/publications/m/item/WHO-UCN-GMP-MPAG-2022.01](http://www.who.int/publications/m/item/WHO-UCN-GMP-MPAG-2022.01)
3. Luzzatto L, Ally M, Notaro R. Glucose-6-phosphate dehydrogenase deficiency. *Blood*. 2020;136(11):1225-1240. doi:10.1182/blood.2019000944
4. Cappellini MD, Fiorelli G. Glucose-6-phosphate dehydrogenase deficiency. *Lancet*. 2008;371:64-67
5. Luzzatto L, Seneca E. G6PD deficiency: a classic example of pharmacogenetics with on-going clinical implications. *Br J Haematol*. 2014;164:469-480
4. OMIM. 305900 Glucose-6-phosphate dehydrogenase. Johns Hopkins University; 1987. Updated April 28, 2023.

---

Accessed July 3, 2023. Available at [www.omim.org/entry/305900](http://www.omim.org/entry/305900)

5. Relling MV, McDonagh EM, Chang T, et al. Clinical Pharmacogenetics Implementation Consortium (CPIC) guidelines for rasburicase therapy in the context of G6PD deficiency genotype. *Clin Pharmacol Ther.* 2014 Aug;96(2):169-174

## Performance

### Method Description

Genomic DNA is extracted from whole blood. The *G6PD* gene is amplified by polymerase chain reaction (PCR). The PCR products are then purified and sequenced in both directions using fluorescent dye-terminator chemistry. Sequencing products are separated on an automated sequencer and trace files analyzed for variations in the exons and intron/exon boundaries of all exons using variant detection software and visual inspection. Variant nomenclature is based on GenBank accession number NM\_001042351.2 using human genome assembly GRCh37 (hg19). (Unpublished Mayo method)

### PDF Report

No

### Day(s) Performed

Monday, Wednesday

### Report Available

3 to 7 days

### Specimen Retention Time

Whole blood/Saliva: 2 weeks; Extracted DNA: 3 months

### Performing Laboratory Location

Rochester

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

### CPT Code Information

81249

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
G6PDZ	G6PD Full Gene Sequencing, V	94231-8

Result ID	Test Result Name	Result LOINC® Value
618837	G6PD Phenotype	47998-0
618838	Result Details	82939-0
618839	Interpretation	69047-9
618840	Additional Information	48767-8
618841	Method	85069-3
618842	Disclaimer	62364-5
618843	Reviewed By	18771-6