

## Overview

### Useful For

Determining the cause for hyper- or hypokalemia

### Special Instructions

- [Urine Preservatives-Collection and Transportation for 24-Hour Urine Specimens](#)

### Method Name

Potentiometric, Indirect Ion-Selective Electrode (ISE)

### NY State Available

Yes

## Specimen

### Specimen Type

Urine

### Necessary Information

**24-Hour volume (in milliliters) is required.**

### Specimen Required

**Supplies:** Sarstedt 5 mL Aliquot Tube (T914)

**Collection Container/Tube:** 24-Hour graduated urine container with no metal cap or glued insert

**Submission Container/Tube:** Plastic, 5 mL tube or a clean, plastic aliquot container with no metal cap or glued insert

**Specimen Volume:** 5 mL

#### Collection Instructions:

1. Collect urine for 24 hours.
2. Refrigerate specimen within 4 hours of completion of 24-Hour collection.

#### Additional Information:

See [Urine Preservatives-Collection and Transportation for 24-Hour Urine Specimens](#) for multiple collections.

### Urine Preservative Collection Options

**Note:** The addition of preservative or application of temperature controls **must occur within 4 hours of completion** of the collection.

|                 |           |
|-----------------|-----------|
| Ambient         | OK        |
| Refrigerate     | Preferred |
| Frozen          | OK        |
| 50% Acetic Acid | OK        |
| Boric Acid      | OK        |

|                      |    |
|----------------------|----|
| Diazolidinyl Urea    | OK |
| 6M Hydrochloric Acid | OK |
| 6M Nitric Acid       | No |
| Sodium Carbonate     | OK |
| Thymol               | OK |
| Toluene              | No |

### Specimen Minimum Volume

1 mL

### Reject Due To

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

### Specimen Stability Information

| Specimen Type | Temperature              | Time    | Special Container |
|---------------|--------------------------|---------|-------------------|
| Urine         | Refrigerated (preferred) | 14 days |                   |
|               | Ambient                  | 7 days  |                   |
|               | Frozen                   | 30 days |                   |

## Clinical & Interpretive

### Clinical Information

Potassium (K<sup>+</sup>) is the major intracellular cation. Functions of potassium include regulation of neuromuscular excitability, heart contractility, intracellular fluid volume, and hydrogen ion concentration. The physiologic function of K<sup>+</sup> requires the body to maintain a low extracellular fluid concentration of the cation; the intracellular concentration is 20 times greater than the extracellular K<sup>+</sup> concentration. Only 2% of total body K<sup>+</sup> circulates in the plasma.

The kidneys provide the most important regulation of K<sup>+</sup>. The proximal tubules reabsorb almost all the filtered K<sup>+</sup>. Under the influence of aldosterone, the remaining K<sup>+</sup> can then be secreted into the urine in exchange for sodium in both the collecting ducts and the distal tubules. Thus, the distal nephron is the principal determinant of urinary K<sup>+</sup> excretion.

Decreased excretion of K<sup>+</sup> in acute kidney disease and end-stage kidney failure are common causes of prolonged hyperkalemia.

Renal losses of K<sup>+</sup> may occur during the diuretic (recovery) phase of acute tubular necrosis, during administration of non-potassium sparing diuretic therapy, and during states of excess mineralocorticoid or glucocorticoid.

### Reference Values

> or =18 years: 16-105 mmol/24 hours

Reference values have not been established for patients who are younger than 18 years of age.

---

**Interpretation**

Hypokalemia reflecting true total body deficits of potassium (K<sup>+</sup>) can be classified into renal and nonrenal losses based on the daily excretion of K<sup>+</sup> in the urine.

During hypokalemia, if urine excretion of K<sup>+</sup> is below 30 mEq/day, it can be concluded that kidney reabsorption of K<sup>+</sup> is appropriate. In this situation, the causes for the hypokalemic state are either decreased K<sup>+</sup> intake or extra renal loss of K<sup>+</sup> rich fluid.

Urine excretion of more than 30 mEq/d in a hypokalemia setting is inappropriate and indicates that the kidneys are the primary source of the lost K<sup>+</sup>.

**Cautions**

Ion-selective electrodes are selective for the ion in question but are not specific. Other monovalent cations may interfere but not in the physiologic range.

**Clinical Reference**

1. Delaney MP, Lamb EJ: Kidney disease. In: Rifai N, Horvath AR, Wittwer CT, eds. Tietz Textbook of Clinical Chemistry and Molecular Diagnostics. 6th ed. Elsevier; 2018:1308-1309
2. Toffaletti J: Electrolytes. In: Dufour DR, Rifai N, eds. Professional Practice in Clinical Chemistry: A Review. AACCC Press; 1993

**Performance****Method Description**

The ion-selective electrode (ISE) module performs indirect measurement of electromotive force (EMF). The ISE module measures the EMF difference between an ion-selective electrode and a reference electrode. The EMF of the ion-selective electrode is dependent on the ion concentration of the sample. The EMF of the reference electrode is constant. An electronic calculation circuit converts EMF of the sample to the ion concentration of the sample. (Package insert: ISE. Roche Diagnostics; V14.0 02/2018)

**PDF Report**

No

**Day(s) Performed**

Monday through Sunday

**Report Available**

Same day/1 to 2 days

**Specimen Retention Time**

7 days

**Performing Laboratory Location**

Mayo Clinic Laboratories - Rochester Main Campus

---

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test has been cleared, approved, or is exempt by the US Food and Drug Administration and is used per manufacturer's instructions. Performance characteristics were verified by Mayo Clinic in a manner consistent with CLIA requirements.

**CPT Code Information**

84133

**LOINC® Information**

| Test ID | Test Order Name     | Order LOINC® Value |
|---------|---------------------|--------------------|
| KUR     | Potassium, 24 HR, U | 2829-0             |

| Result ID | Test Result Name    | Result LOINC® Value |
|-----------|---------------------|---------------------|
| KURU      | Potassium, 24 Hr, U | 2829-0              |
| TM28      | Collection Duration | 13362-9             |
| VL26      | Urine Volume        | 3167-4              |