

Hypertriglyceridemia Gene Panel, Varies

Overview

Useful For

Providing a genetic evaluation for patients with a personal or family history suggestive of hereditary forms of primary hypertriglyceridemia and related conditions.

Genetics Test Information

This test utilizes next-generation sequencing to detect single nucleotide and copy number variants in 13 genes associated with primary hypertriglyceridemia and related conditions: *APOA5, APOC2, APOE, CREB3L3, GPD1, GPIHBP1, LCAT, LIPA, LIPC, LMF1, LPL,* and *LRP6*. See <u>Targeted Genes and Methodology Details for Hypertriglyceridemia Gene</u> Panel and Method Description for additional details.

Identification of a disease-causing variant may assist with diagnosis, prognosis, clinical management, familial screening, and genetic counseling for hypertriglyceridemia and related conditions.

This test also reports homozygous status for the *APOE* E2 allele, a risk allele for type III hyperlipoproteinemia. This test does NOT report other, non-cardiovascular *APOE* disease associations, including Alzheimer disease.

Prior Authorization is available for this assay.

Special Instructions

- Informed Consent for Genetic Testing
- Informed Consent for Genetic Testing (Spanish)
- Hereditary Dyslipidemia Patient Information
- Targeted Genes and Methodology Details for Hypertriglyceridemia Gene Panel
- Hypertriglyceridemia Gene Panel (HYPTG) Prior Authorization Ordering Instructions

Method Name

Sequence Capture and Targeted Next-Generation Sequencing followed by Polymerase Chain Reaction (PCR) and Sanger Sequencing

NY State Available

Yes

Specimen

Specimen Type

Varies

Ordering Guidance

Customization of this panel and single gene analysis for any gene present on this panel are available. or more



Hypertriglyceridemia Gene Panel, Varies

information, see CGPH / Custom Gene Panel, Hereditary, Next-Generation Sequencing, Varies.

Targeted testing for familial variants (also called site-specific or known mutations testing) is available for the genes on this panel. See FMTT / Familial Variant, Targeted Testing, Varies. To obtain more information about this testing option, call 800-533-1710.

Shipping Instructions

Specimen preferred to arrive within 96 hours of collection.

Necessary Information

<u>Prior Authorization</u> is available, **but not required**, for this test. If proceeding with the prior authorization process, submit the required form with the specimen.

Specimen Required

Patient Preparation: A previous bone marrow transplant from an allogenic donor will interfere with testing. Call

800-533-1710 for instructions for testing patients who have received a bone marrow transplant.

Specimen Type: Whole blood

Container/Tube:

Preferred: Lavender top (EDTA) or yellow top (ACD)

Acceptable: Any anticoagulant Specimen Volume: 3 mL Collection Instructions:

1. Invert several times to mix blood.

2. Send whole blood specimen in original tube. Do not aliquot.

Specimen Stability Information: Ambient (preferred)/Refrigerated

Forms

1. **New York Clients-Informed consent is required.** Document on the request form or electronic order that a copy is on file

The following documents are available:

- -Informed Consent for Genetic Testing (T576)
- -Informed Consent for Genetic Testing (Spanish) (T826)
- 2. Hereditary Dyslipidemia Patient Information
- 3. Hypertriglyceridemia Gene Panel (HYPTG) Prior Authorization Ordering Instructions
- 4. If not ordering electronically, complete, print, and send a <u>Cardiovascular Test Request</u> (T724) with the specimen.

Specimen Minimum Volume

1 mL

Reject Due To

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|-------------|------|-------------------|
| Varies | Varies | | |



Hypertriglyceridemia Gene Panel, Varies

Clinical & Interpretive

Clinical Information

Hypertriglyceridemia (HTG), or abnormally elevated triglyceride concentration in the blood, is present in approximately 30% of adults in the United States.(1) HTG is frequently associated with other abnormalities such as abdominal obesity, insulin resistance, low high-density lipoprotein (HDL), and hypertension, which are linked to coronary artery disease and metabolic syndrome. Severe hypertriglyceridemia is associated with an increased risk of acute pancreatitis.

Hypertriglyceridemia can be classified into primary and secondary types. Primary hypertriglyceridemia accounts for less than 5% of cases and is due to rare, monogenic conditions with disease-causing variants resulting in disordered triglyceride metabolism.(2) The majority of hypertriglyceridemia is secondary hypertriglyceridemia, due to a combination of lifestyle factors such as high fat diet, obesity, diabetes, hypothyroidism, or certain medications, in addition to various common, low impact genetic variants that cumulatively have an impact on triglyceride metabolism.

This test includes analysis of several genes associated with monogenic forms of primary hypertriglyceridemia and related conditions. Autosomal recessive conditions tested for by this panel include combined lipase deficiency (*LMF1*), hepatic lipase deficiency (*LIPC*), lysosomal acid lipase deficiency (also known as Wolman disease) and cholesteryl ester storage disease (*LIPA*), lecithin:cholesterol acyltransferase deficiency (also known as Norum disease or fish-eye disease) (*LCAT*), hyperlipoproteinemia type 1b (*APOC2*), transient infantile hypertriglyceridemia (*GPD1*), hyperlipoproteinemia type ID (*GPIHBP1*), and lipoprotein lipase deficiency (*LPL*). Autosomal dominant conditions tested for by this panel include familial hypertriglyceridemia and late onset hyperchylomicronemia (*APOA5*), hypertriglyceridemia 2 (*CREB3L3*), and familial combined hyperlipidemia (*LPL*). In addition, this test reports homozygous status for the *APOE* E2 allele, a risk allele for type III hyperlipoproteinemia.

Reference Values

An interpretive report will be provided.

Interpretation

All detected variants are evaluated according to American College of Medical Genetics and Genomics recommendations.(3) Variants are classified based on known, predicted, or possible pathogenicity and reported with interpretive comments detailing their potential or known significance.

Cautions

Clinical Correlations:

Test results should be interpreted in the context of clinical findings, family history, and other laboratory data. Misinterpretation of results may occur if the information provided is inaccurate or incomplete.

If testing was performed because of a clinically significant family history, it is often useful to first test an affected family member. Detection of a reportable variant in an affected family member would allow for more informative testing of at-risk individuals.

To discuss the availability of additional testing options or for assistance in the interpretation of these results, contact the Mayo Clinic Laboratories genetic counselors at 800-533-1710.



Hypertriglyceridemia Gene Panel, Varies

Technical Limitations:

Next-generation sequencing may not detect all types of genomic variants. In rare cases, false-negative or false-positive results may occur. The depth of coverage may be variable for some target regions; assay performance below the minimum acceptable criteria or for failed regions will be noted. Given these limitations, negative results do not rule out the diagnosis of a genetic disorder. If a specific clinical disorder is suspected, evaluation by alternative methods can be considered.

There may be regions of genes that cannot be effectively evaluated by sequencing or deletion and duplication analysis as a result of technical limitations of the assay, including regions of homology, high guanine-cytosine (GC) content, and repetitive sequences. Confirmation of select reportable variants will be performed by alternate methodologies based on internal laboratory criteria.

This test is validated to detect 95% of deletions up to 75 base pairs (bp) and insertions up to 47 bp. Deletions-insertions (delins) of 40 or more bp, including mobile element insertions, may be less reliably detected than smaller delins.

Deletion/Duplication Analysis:

This analysis targets single and multi-exon deletions/duplications; however, in some instances single exon resolution cannot be achieved due to isolated reduction in sequence coverage or inherent genomic complexity. Balanced structural rearrangements (such as translocations and inversions) may not be detected.

This test is not designed to detect low levels of mosaicism or to differentiate between somatic and germline variants. If there is a possibility that any detected variant is somatic, additional testing may be necessary to clarify the significance of results.

Genes may be added or removed based on updated clinical relevance. For the most up to date list of genes included in this test and detailed information regarding gene specific performance and technical limitations, see Method Description or contact a laboratory genetic counselor.

If the patient has had an allogeneic hematopoietic stem cell transplant or a recent blood transfusion, results may be inaccurate due to the presence of donor DNA. Call Mayo Clinic Laboratories for instructions for testing patients who have received a bone marrow transplant.

Reclassification of Variants:

At this time, it is not standard practice for the laboratory to systematically review previously classified variants on a regular basis. The laboratory encourages healthcare providers to contact the laboratory at any time to learn how the classification of a particular variant may have changed over time.

Variant Evaluation:

Evaluation and categorization of variants are performed using published American College of Medical Genetics and Genomics and the Association for Molecular Pathology recommendations as a guideline. (3) Other gene-specific guidelines may also be considered. Variants are classified based on known, predicted, or possible pathogenicity and reported with interpretive comments detailing their potential or known significance. Variants classified as benign or likely benign are not reported.

Multiple in silico evaluation tools may be used to assist in the interpretation of these results. The accuracy of predictions



Hypertriglyceridemia Gene Panel, Varies

made by in silico evaluation tools is highly dependent upon the data available for a given gene, and periodic updates to these tools may cause predictions to change over time. Results from in silico evaluation tools should be interpreted with caution and professional clinical judgment.

Rarely, incidental or secondary findings may implicate another predisposition or presence of active disease. Incidental findings may include, but are not limited to, results related to the sex chromosomes. These findings will be carefully reviewed to determine whether they will be reported.

Clinical Reference

- 1. Pejic RN and Lee DT. Hypertriglyceridemia. J Am Board Fam Med. 2006;19(3):310-316
- 2. Yuan G, Al-Shali KZ, Hegele RA. Hypertriglyceridemia: its etiology, effects and treatment. CMAJ. 2007;176(8): 1113-1120
- 3. Richards S, Aziz N, Bale S, et al. Standards and guidelines for the interpretation of sequence variants: a joint consensus recommendation of the American College of Medical Genetics and Genomics and the Association for Molecular Pathology. Genet Med. 2015;17(5):405-424

Performance

Method Description

Next-generation sequencing (NGS) and/or Sanger sequencing are performed to test for the presence of variants in coding regions and intron/exon boundaries of the genes analyzed, as well as some other regions that have known disease-causing variants. The human genome reference GRCh37/hg19 build was used for sequence read alignment. At least 99% of the bases are covered at a read depth over 30X. Sensitivity is estimated at above 99% for single nucleotide variants, above 94% for deletions-insertions (delins) less than 40 base pairs (bp), above 95% for deletions up to 75 bp and insertions up to 47 bp. NGS and/or a polymerase chain reaction-based quantitative method is performed to test for the presence of deletions and duplications in the genes analyzed.

There may be regions of genes that cannot be effectively evaluated by sequencing or deletion and duplication analysis as a result of technical limitations of the assay, including regions of homology, high guanine-cytosine (GC) content, and repetitive sequences. See <u>Targeted Genes and Methodology Details for Hypertriglyceridemia Gene Panel</u> for details regarding the targeted genes analyzed for each test and specific gene regions not routinely covered.(Unpublished Mayo method)

Confirmation of select reportable variants may be performed by alternate methodologies based on internal laboratory criteria.

Genes analyzed: APOA5, APOC2, APOE, CREB3L3, GPD1, GPIHBP1, LCAT, LIPA, LIPC, LMF1, LPL, and LRP6

PDF Report

Supplemental

Day(s) Performed

Varies



Hypertriglyceridemia Gene Panel, Varies

Report Available

28 to 42 days

Specimen Retention Time

Whole blood: 2 weeks (if available); Extracted DNA: 3 months

Performing Laboratory Location

Mayo Clinic Laboratories - Rochester Main Campus

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

CPT Code Information

81479

Prior Auhtorization

Insurance preauthorization is available for this testing; forms are available.

Patient financial assistance may be available to those who qualify. Patients who receive a bill from Mayo Clinic Laboratories will receive information on eligibility and how to apply.

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|---------------------------------|--------------------|
| HYPTG | Hypertriglyceridemia Gene Panel | 51966-0 |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------------|---------------------|
| 617324 | Test Description | 62364-5 |
| 617325 | Specimen | 31208-2 |
| 617326 | Source | 31208-2 |
| 617327 | Result Summary | 50397-9 |
| 617328 | Result | 82939-0 |
| 617329 | Interpretation | 69047-9 |
| 617330 | Additional Results | 82939-0 |
| 617331 | Resources | 99622-3 |
| 617332 | Additional Information | 48767-8 |



Hypertriglyceridemia Gene Panel, Varies

| 617333 | Method | 85069-3 |
|--------|----------------|---------|
| 617334 | Genes Analyzed | 48018-6 |
| 617335 | Disclaimer | 62364-5 |
| 617336 | Released By | 18771-6 |