

Hematologic Disorders, DNA Extract and Hold,
Varies

#### Overview

#### **Useful For**

Reserving DNA on any specimen for which molecular analysis requiring DNA may be necessary at a future date, ensuring that adequate material for testing is available

#### **Testing Algorithm**

DNA will be extracted from the specimen and stored at -80 degrees Celsius for one year from the date of the extraction process was performed. A report of "Performed" will be sent and a \$50.00 processing fee will be assessed. To request molecular testing from the Molecular Hematopathology Laboratory's test menu, call 800-533-1710. Additional charges will be assessed for testing performed.

# **Special Instructions**

• Hematopathology Patient Information

# **Method Name**

**Nucleic Acid Extraction** 

#### **NY State Available**

Yes

# Specimen

## **Specimen Type**

Varies

#### Specimen Required

Specimen must arrive within 7 days of collection. Collect and package specimen as close to shipping time as possible.

Submit only 1 of the following specimens:

**Specimen Type**: Blood **Container/Tube**:

Preferred: Lavender top (EDTA)
Acceptable: Yellow top (ACD)
Specimen Volume: 4 mL
Collection Instructions:

- 1. Invert several times to mix blood.
- 2. Send whole blood specimen in original tube. Do not aliquot.
- 3. Label specimen as blood.



Hematologic Disorders, DNA Extract and Hold,
Varies

Specimen Type: Bone marrow

Container/Tube:

Preferred: Lavender top (EDTA)
Acceptable: Yellow top (ACD)
Specimen Volume: 2 mL
Collection Instructions:

- 1. Invert several times to mix bone marrow.
- 2. Send bone marrow specimen in original tube. **Do not aliquot**.
- 3. Label specimen as bone marrow.

#### **Forms**

- 1. Hematopathology Patient Information (T676)
- 2. If not ordering electronically, complete, print, and send a <u>Hematopathology/Cytogenetics Test Request</u> (T726) with the specimen.

# **Specimen Minimum Volume**

Blood, Bone Marrow: 1 mL

## **Reject Due To**

Gross	Reject
hemolysis	
Paraffin-embe	Reject
dded material	
Fluids	
Heparinized	
samples	

# **Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)	7 days	PURPLE OR PINK TOP/EDTA
	Refrigerated	7 days	PURPLE OR PINK TOP/EDTA

# **Clinical & Interpretive**

#### **Clinical Information**

It is frequently useful to obtain DNA from clinical samples containing a hematopoietic neoplasm at the time of diagnosis so appropriate material is available for molecular analysis should subsequent testing be required. For example, when a diagnosis of acute myelogenous leukemia is made, there is a delay before karyotype information, which determines whether testing for molecular prognostic markers is necessary, is available. After this delay, the diagnostic sample is usually no longer available, or the DNA has degraded to such an extent that it is no longer adequate for testing. Thus, it



Hematologic Disorders, DNA Extract and Hold,
Varies

is useful to obtain DNA on such specimens promptly at diagnosis and retain it until it is known whether additional testing is necessary.

### **Reference Values**

Not applicable

# Interpretation

A report of "Performed" will be sent and a \$50 processing fee will be assessed. No interpretation will be given. Should the sample be used in future testing, interpretation would be incorporated with the final testing.

#### **Cautions**

DNA will be held at -80 degrees Celsius for one year from the time of extraction. The DNA may only be used for suitable tests offered by the Molecular Hematopathology Laboratory at Mayo Clinic. It cannot be shipped to another institution for testing and cannot be released for research use.

Only liquid blood specimens will be processed. Paraffin-embedded material and fresh, solid tissue will not be accepted at this time.

#### **Performance**

### **Method Description**

DNA is extracted using an automated system. The system includes a sample lysis step followed by binding of nucleic acid to a solid phase for separation from other cellular components and a final nucleic acid elution step.

#### PDF Report

No

#### Day(s) Performed

Monday through Saturday

# **Report Available**

5 days

#### **Specimen Retention Time**

Blood/Bone marrow: 2 weeks; Extracted DNA:1 year

# **Performing Laboratory Location**

Mayo Clinic Laboratories - Rochester Main Campus

## **Fees & Codes**



Hematologic Disorders, DNA Extract and Hold,
Varies

#### **Fees**

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

### **Test Classification**

Not Applicable

# **LOINC®** Information

Test ID	Test Order Name	Order LOINC® Value
EXHD	DNA Extract and Hold	No LOINC Needed

Result ID	Test Result Name	Result LOINC® Value
MP034	Specimen Type	31208-2
41995	DNA Extract and Hold Result	No LOINC Needed
38305	DNA Extraction	No LOINC Needed