

## Overview

### Useful For

Recurrent common chromosome abnormalities in patients with small lymphocytic lymphoma (SLL)

Distinguishing patients with 11;14 translocations who have mantle cell lymphoma (MCL) from patients who have SLL

Detecting patients with atypical SLL or other forms of B-cell lymphoma associated with translocations between *IGH* and *BCL3*

### Reflex Tests

Test Id	Reporting Name	Available Separately	Always Performed
_PRAG	Probe, Each Additional (SLL)	No, (Bill Only)	No

### Testing Algorithm

This test includes a charge for the probe application, analysis, and professional interpretation of results for one probe set (2 individual fluorescence in situ hybridization [FISH] probes). Additional charges will be incurred for all reflex or additional probes performed. No analysis charges will be incurred if an insufficient number of representative cells are available for analysis.

This test may be ordered in 2 distinct ways allowing different combinations of probes to be analyzed based on the clinical question, including the standard small lymphocytic lymphoma (SLL) FISH panel and the individual SLL FISH probes (per client request).

**If the patient is being evaluated for known abnormalities, targeted probes must be listed in the probe request field.** If no specific panel or FISH probes are indicated, the standard panel will be performed.

The standard SLL FISH panel includes testing for the following abnormalities, using the FISH probes listed:

- 6q-, D6Z1/MYB
- 11q-, D11Z1/ATM
- +12, D12Z3/MDM2
- 13q-, D13S319/LAMP1
- 17p-, TP53/D17Z1
- t(11;14), CCND1::IGH

When 3 IGH signals are identified suggesting an *IGH* rearrangement and no fusion with CCND1 is observed, additional testing with the t(14;19)(q32;q13) IGH::BCL3 FISH probe will be performed.

Appropriate ancillary probes may be performed at consultant discretion to render comprehensive assessment. Any additional probes will have the results included within the final report and will be performed at an additional charge.

---

**Method Name**

Fluorescence In Situ Hybridization (FISH)

**NY State Available**

Yes

**Specimen****Specimen Type**

Tissue

**Ordering Guidance**

This test does not include a pathology consultation. If a pathology consultation is requested, order PATHC / Pathology Consultation and the appropriate fluorescence in situ hybridization test (FISH) test will be added and performed at an additional charge.

Mayo Hematopathology Consultants are involved in both the preanalytic (tissue adequacy and probe selection, when applicable) and postanalytic (interpretation of FISH results in context of specific case, when applicable) phases.

This test is **not appropriate** for testing blood and bone marrow from patients with chronic lymphocytic leukemia. See CLLDF / Chronic Lymphocytic Leukemia (CLL), Diagnostic FISH, Varies or CLLMF / Chronic Lymphocytic Leukemia (CLL), Specified FISH, Varies.

**Shipping Instructions**

Advise Express Mail or equivalent if not on courier service.

**Necessary Information**

**1. A pathology report is required for testing to be performed.** If not provided, appropriate testing and/or interpretation may be compromised or delayed. Acceptable pathology reports include working drafts, preliminary pathology, or surgical pathology reports.

**2. The following information must be included in the report provided:**

- Patient name
- Block number- must be on all blocks, slides, and paperwork
- Date of collection
- Tissue source

**3. A reason for testing must be provided.** If this information is not provided, an appropriate indication for testing may be entered by Mayo Clinic Laboratories.

**Specimen Required**

**Submit only 1 of the following specimens:**

**Preferred**

**Specimen Type:** Tissue block

**Collection Instructions:** Submit a formalin-fixed, paraffin-embedded tumor tissue block. Blocks prepared with alternative fixation methods may be acceptable; provide fixation method used.

**Acceptable**

**Specimen Type:** Tissue slides

**Slides:** 1 Hematoxylin and eosin-stained (H and E) stained and 10 unstained

**Collection Instructions:** Submit 1 slide stained with H and E and 10 consecutive, unstained, positively charged, unbaked slides with 5-micron-thick sections of the tumor tissue.

## Forms

If not ordering electronically, complete, print, and send a [Hematopathology/Cytogenetics Test Request](#) (T726) with the specimen.

## Specimen Minimum Volume

Slides: 1 Hematoxylin and eosin-stained and 6 unstained

## Reject Due To

All specimens will be evaluated at Mayo Clinic Laboratories for test suitability.

## Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Tissue	Ambient (preferred)		
	Refrigerated		

## Clinical & Interpretive

### Clinical Information

Small lymphocytic lymphoma (SLL) is the nonleukemic form of chronic lymphocytic leukemia (CLL), one of the most common leukemias in adults. The most frequently seen cytogenetic abnormalities in SLL involve chromosomes 6, 11, 12, 13 and 17. These are detected and quantified using the SLL fluorescence in situ hybridization (FISH) panel.

Cytogenetics has proven to be a reliable predictor of outcome for patients with CLL. It is unknown if SLL has the same prognostic significance when these genetic abnormalities are observed.

This FISH test detects an abnormal clone in approximately 65% of patients with SLL. Patients with t(11;14)(q13;q32) associated with *CCND1::IGH* fusion, have mantle cell lymphoma which can be distinguished from SLL and other B-cell lymphomas with this assay. Patients with t(14;19)(q32;q13.3) associated with *IGH::BCL3* fusion, may have an atypical form of SLL or another B-cell lymphoma.

### Reference Values

An interpretive report will be provided.

### Interpretation

A neoplastic clone is detected when the percent of cells with an abnormality exceeds the normal reference range for any given probe set.

A positive result is not diagnostic for small lymphocytic lymphoma but may provide relevant prognostic information.

The absence of an abnormal clone does not rule out the presence of an SLL clone or another neoplastic disorder.

**Cautions**

This test is not approved by the US Food and Drug Administration, and it is best used as an adjunct to existing clinical and pathologic information.

Fixatives other than formalin (eg, Prefer, Bouin's) may not be successful for fluorescence in situ hybridization (FISH) assays. Non-formalin fixed samples will not be rejected.

Paraffin-embedded tissues that have been decalcified may not be successful for FISH analysis. The success rate of FISH studies on decalcified tissue is approximately 50%, but FISH will be attempted if sufficient tumor is present for analysis.

If no FISH signals are observed post hybridization, the case will be released indicating a lack of FISH results.

**Clinical Reference**

1. Swerdlow SH, Campo E, Harris NL eds. WHO Classification of Tumours of Haematopoietic and Lymphoid Tissues. IARC; 2017. WHO Classification of Tumours, Vol 2
2. Shanafelt TD. Predicting clinical outcome in CLL: how and why. Hematology Am Soc Hematol Educ Program. 2009;421-429
3. Van Dyke DL, Werner L, Rassenti LZ, et al. The Dohner fluorescence in situ hybridization prognostic classification of chronic lymphocytic leukaemia (CLL): the CLL Research Consortium experience. Br J Haematol. 2016;173(1):105-113

**Performance****Method Description**

This test is performed using commercially available and laboratory-developed probes. Deletion of chromosomes 6q, 11q, 13q, and 17p, and trisomy of chromosome 12 are detected using enumeration strategy probes. A dual-color, dual-fusion (D-FISH) strategy probe set is used to detect *CCND1::IGH* rearrangements and, for reflex testing, to identify *IGH::BCL3* rearrangements. Paraffin-embedded tissues are cut at 5 microns and mounted on positively charged glass slides. The selection of tissue and the identification of target areas on the hematoxylin and eosin (H and E)-stained slide is performed by a pathologist. Using the H and E-stained slide as a reference, target areas are etched with a diamond-tipped engraving tool on the back of the unstained slide to be assayed. For each probe set, the probes are hybridized to the appropriate target areas and 2 technologists each analyze 50 interphase nuclei (100 total). All results are expressed as the percent abnormal nuclei. (Unpublished Mayo method)

**PDF Report**

No

**Day(s) Performed**

[Monday through Friday](#)

**Report Available**

7 to 10 days

**Specimen Retention Time**

Slides and H and E used for analysis are retained by the laboratory in accordance to with CAP and NYS regulatory requirements. Client provided paraffin blocks and extra unstained slides (if provided) will be returned after testing is complete.

### Performing Laboratory Location

Rochester

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

### CPT Code Information

- 88377-if 1 probe set
- 88377 x 2-if 2 probe sets
- 88377 x 3-if 3 probe sets
- 88377 x 4-if 4 probe sets
- 88377 x 5-if 5 probe sets
- 88377 x 6-if 6 probe sets
- 88377 x 7-if 7 probe sets
- 88377 x 8-if 8 probe sets

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
SLL	SLL, FISH, Tissue	103621-9

Result ID	Test Result Name	Result LOINC® Value
603129	Result Summary	50397-9
603130	Interpretation	69965-2
603131	Result Table	93356-4
603132	Result	62356-1
GC038	Reason for Referral	42349-1
603133	Specimen	31208-2
603134	Source	31208-2
603135	Tissue ID	80398-1
603136	Method	85069-3
603137	Additional Information	48767-8

---

603138	Disclaimer	62364-5
603139	Released By	18771-6