

## Overview

### Method Name

Quantitative Spectrophotometry

### NY State Available

Yes

## Specimen

### Specimen Type

Peritoneal

### Specimen Required

**Specimen Type:** Peritoneal fluid (Ascites, Paracentesis)

**Container/Tube:** Standard transport container

**Specimen volume:** 0.5 mL

**Collection Instructions:** Collect Peritoneal Fluid in a leak-proof container. Centrifuge specimen at room temperature, transfer 0.5 mL peritoneal fluid to plastic vial and Ship frozen. The specimen must remain frozen until received at the performing lab.

### Specimen Minimum Volume

0.2 mL

### Reject Due To

Thawing**	Cold OK; Warm reject
Other/Tissue/Swab	Whole blood, Bronchoalveolar lavage (BAL) specimens, Turbid specimen

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Peritoneal	Refrigerated	7 days	
	Frozen (preferred)	30 days	

## Clinical & Interpretive

### Reference Values

0 - 30 U/L

**Performance****PDF Report**

No

**Day(s) Performed**

Sunday, Tuesday, Thursday

**Report Available**

1 to 8 days

**Performing Laboratory Location**

ARUP Laboratories

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the U.S. Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

**CPT Code Information**

84311

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FADFL	Adenosine Deaminase Peritoneal Fld	49759-4

Result ID	Test Result Name	Result LOINC® Value
FADFL	Adenosine Deaminase Peritoneal Fld	49759-4