

Overview

Reflex Tests

Test Id	Reporting Name	Available Separately	Always Performed
FA13R	Factor 13 1:1 Mix	No	No

Testing Algorithm

If the factor 13 qualitative result states "Lysis", then factor 13 1:1 mix will be performed at an additional charge.

Method Name

Qualitative Solubility Assay

NY State Available

Yes

Specimen

Specimen Type

Plasma Na Cit

Specimen Required

**Collection Container/Tube:** Light-blue top (3.2% sodium citrate)

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 2 mL Platelet-poor plasma

**Collection Instructions:**

- 1. Centrifuge immediately and aliquot 2 mL platelet-poor plasma in a plastic vial.
- 2. Send frozen.

**Additional Information:** Critical frozen. Separate specimens must be submitted when multiple tests are ordered.

Specimen Minimum Volume

Platelet-poor plasma: 1 mL

Reject Due To

Hemolysis	Reject
Clotted specimens	Reject
Serum	Reject
EDTA plasma	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Plasma Na Cit	Frozen	14 days	

Clinical & Interpretive

Clinical Information

Refer to <https://ltd.aruplab.com/>.

Reference Values

Factor XIII, Qualitative: No Lysis

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

4 to 5 days

Performing Laboratory Location

ARUP Laboratories

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

85291-Factor 13, Qualitative  
85291-Factor 13, 1:1 Mix (if appropriate)

# Test Definition: FA13Q

Factor 13, Qualitative, with Reflex to Factor 13  
1:1 Mix

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FA13Q	Factor 13 Qual, Reflex to 1:1 Mix	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z6312	Factor XIII, Qualitative	3241-7