

Overview

Reflex Tests

| Test Id | Reporting Name | Available Separately | Always Performed |
|---------|----------------------------------|----------------------|------------------|
| FSGCU | Specific Gravity Confirmation, U | No | No |

Testing Algorithm

If Creatinine result is less than 300 or greater than 3000, Specific Gravity Confirmation, Urine (FSGCU) will be performed, if appropriate, at no additional charge.

Method Name

Colorimetry (C)

High Performance Liquid Chromatography/Tandem Mass Spectrometry (LS-MS/MS)

Refractometer (REF), if appropriate

NY State Available

Yes

Specimen

Specimen Type

Urine

Specimen Required

Must submit one specimen per order. Specimens cannot be shared between multiple orders.

Specimen Type: Urine

Container/Tube: Plastic preservative-free urine container

Specimen Volume: 4 mL

Collection Instructions: Collect 4 mL random urine without preservative. Send specimen refrigerated in a plastic (preservative-free) urine container.

Specimen Minimum Volume

1.6 mL

Reject Due To

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|--------------------------|---------|-------------------|
| Urine | Refrigerated (preferred) | 30 days | |
| | Frozen | 14 days | |

Clinical & Interpretive

Reference Values

Reporting limit determined each analysis

Creatinine (mg/L)

U.S. Population (10th-90th percentiles, median)

All participants:

335 - 2370 mg/L, median: 1180 (n=22,245)

Males:

495 - 2540 mg/L, median: 1370 (n=10, 610)

Females:

273 - 2170 mg/L, median: 994 (n=11,635)

Thiosulfate (mcg/mL)

Generally less than 9.2 mcg/mL (based on a median creatinine concentration of 1.18 g/L)

Thiosulfate (Creatinine corrected) (mcg/g Creat)

Generally less than 7.8 mg/g creatinine

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

9 to 16 days

Specimen Retention Time

2 weeks

Performing Laboratory Location

NMS Labs

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

82542

82570

81002 (if appropriate)

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|--------------------|--------------------|
| FSULU | Thiosulfate, Urine | Not Provided |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------------------------|---------------------|
| Z6354 | Creatinine | 2161-8 |
| Z6355 | Thiosulfate | 12420-6 |
| Z6356 | Thiosulfate (Creatinine corrected) | 32274-3 |