

# **Test Definition: HIPA**

Hippuric Acid, Urine

## Overview

#### Method Name

Colorimetry (C) High Performance Liquid Chromatography/Tandem Mass Spectrometry (LC-MS/MS)

#### NY State Available

Yes

Specimen

## Specimen Type

Urine

## Specimen Required

\*\*\*Must submit one specimen per order. Specimens cannot be shared between multiple orders.\*\*\*
Specimen Type: Urine
Container/Tube: Plastic preservative-free urine container
Specimen Volume: 3 mL
Collection Instructions: Collect 3 mL random urine specimen at end of shift. Send specimen refrigerated in a plastic, preservative-free urine container.

## **Specimen Minimum Volume**

1.1 mL

## **Reject Due To**

## Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Urine	Refrigerated (preferred)	30 days	
	Ambient	5 days	
	Frozen	30 days	

## **Clinical & Interpretive**

## **Reference Values**

Reporting limit determined each analysis



# **Test Definition: HIPA**

Hippuric Acid, Urine

Creatinine (mg/L): U.S. Population (10th – 90th percentiles, median) All participants: 335 - 2370 mg/L, median: 1180 (n=22,245) Males: 495 - 2540 mg/L, median: 1370 (n=10,610) Females: 273 - 2170 mg/L, median 994 (n=11,635) Hippuric Acid (g/L) Normal for unexposed populations is generally less than 1.6 g/L.

Hippuric Acid (Creatinine corrected) (g/g Creat) Normal for unexposed populations is generally less than 1.5 g/g creatinine.

Specific Gravity Confirmation Physiologic range: 1.010 - 1.030

## Performance

PDF Report No

Day(s) Performed Monday through Sunday

**Report Available** 9 - 16 days

Specimen Retention Time 2 weeks

Performing Laboratory Location NMS Labs

# Fees & Codes

## Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.



## **Test Classification**

This test was developed and its performance characteristics determined by NMS Labs. It has not been cleared or approved by the U.S. Food and Drug Administration.

## **CPT Code Information**

82570 83921 81002 (if appropriate)

## LOINC<sup>®</sup> Information

Test ID	Test Order Name	Order LOINC <sup>®</sup> Value
НІРА	Hippuric Acid, U Not Provided	
Result ID	Test Result Name	Result LOINC <sup>®</sup> Value

21305	Creatinine	2161-8
Z1341	Reporting Limit	19147-8
Z1309	Hippuric Acid	6709-0
Z1342	Reporting Limit	19147-8
Z1306	Hippuric Acid(Creatinine corrected)	13754-7
Z1343	Reporting Limit	19147-8
Z2513	Specific Gravity Confirmation	In Process