

## Overview

### Useful For

Determining the subspecies of a *Mycobacterium abscessus* complex culture isolate

### Testing Algorithm

When this test is ordered, subspecies identification will always be performed at an additional charge if the organism is not already identified to the subspecies level.

### Special Instructions

- [Infectious Specimen Shipping Guidelines](#)

### Method Name

Next-Generation Sequencing (NGS)

### NY State Available

Yes

## Specimen

### Specimen Type

Varies

### Ordering Guidance

This test is used to identify the subspecies of *Mycobacterium abscessus* from a known *M abscessus* isolate.

For the identification of *M abscessus* (not the subspecies) from culture isolate growth, order CTBID / Culture Referred for Identification, *Mycobacterium* and *Nocardia*, Varies.

### Shipping Instructions

1. For shipping information see [Infectious Specimen Shipping Guidelines](#).
2. Place specimen in a large infectious container and label as an etiologic agent/infectious substance.

### Necessary Information

**Specimen source and suspected organism identification are required.**

### Specimen Required

**Supplies:** Infectious Container, Large (T146)

**Specimen Type:** *Mycobacterium abscessus* isolate growing in pure culture

**Container/Tube:**

**Preferred:** Solid slant medium (eg. Middlebrook 7H10, 7H11 or Lowenstein Jensen agar)

**Acceptable:** Broth medium (eg. Mycobacteria Growth Indicator Tube, 7H9 broth, BACT/ALERT MP or VersaTREK bottle).

**Note:** Broth specimens will require subculture which will delay results.

**Specimen Volume:**

Solid slant medium: Isolate with visible growth

Broth medium: Greater than or equal to 3 mL of broth culture

**Collection Instructions:**

1. Organism must be in pure culture, actively growing. **Do not submit mixed cultures.**
2. Place specimen in a large infectious container and label as an etiologic agent/infectious substance.

**Specimen Minimum Volume**

See Specimen Required

**Reject Due To**

Agar plate	Reject
Mixed culture	Reject

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Varies	Ambient (preferred)		
	Refrigerated		

**Clinical & Interpretive**

**Clinical Information**

This assay provides a subspecies-level identification of microbiologic culture isolates previously identified to be *Mycobacterium abscessus*. Subspecies level identification can be important for patient care and for epidemiologic investigations since the subspecies can differ in macrolide susceptibility.

**Reference Values**

Not applicable

**Interpretation**

This assay can differentiate the subspecies of *Mycobacterium abscessus*, which are, *M abscessus* subspecies *abscessus*, *M abscessus* subspecies *bolletii*, and *M abscessus* subspecies *massiliense*.

**Cautions**

Only pure culture isolates of *Mycobacterium abscessus* should be submitted. Mixed cultures will result in a delay because the *M abscessus* must be isolated prior to performing the next-generation sequencing assay.

This assay has not been verified for the direct detection and direct subspecies identification of *M abscessus* from clinical

specimens. It is intended for use on microbiologic culture isolates already identified as *M abscessus*.

**Clinical Reference**

- 1 Daley CL, Iaccarino JM, Lange C, et al. Treatment of Nontuberculous Mycobacterial Pulmonary Disease: An Official ATS/ERS/ESCMID/IDSA Clinical Practice Guideline. *Clin Infect Dis*. 2020;71(4):905-913. doi:10.1093/cid/ciaa1125
2. Kumar K, Daley CL, Griffith DE, Loebinger MR. Management of *Mycobacterium avium* complex and *Mycobacterium abscessus* pulmonary disease: therapeutic advances and emerging treatments. *Eur Respir Rev*. 2022;31(163):210212. Published 2022 Feb 9. doi:10.1183/16000617.0212-2021

**Performance****Method Description**

The subspecies of *Mycobacterium abscessus* are differentiated using the Genoscreen Deeplex Myc-TB assay. This test is a targeted next generation sequencing test for identification of *Mycobacterium* species using deep sequencing of a single 24-plexed amplicon mix for identification. The gene target used in this assay for species identification is primarily *hsp65* but additional targets such as *rrl*, *rrs*, *rpoB* are used to help with nontuberculous mycobacteria identification to permit discrimination of closely related species and subspecies. This test was validated to demonstrate its ability to differentiate the 3 subspecies within the *M abscessus* from microbiologic culture isolates using the Nextera XT DNA Library Preparation kit, the MiSeq Reagent Kit V2. 300 cycle kit and the Illumina MiSeq sequencing platform. (Buckwalter SP, Olson SL, Fida M, et al. *Mycobacterium abscessus* subspecies identification using the Deeplex Myc-TB targeted NGS assay. *J Clin Microbiol*. 2023;61(10):e0048923. doi:10.1128/jcm.00489-23)

**PDF Report**

No

**Day(s) Performed**

Monday through Friday

**Report Available**

7 to 14 days

**Specimen Retention Time**

1 year

**Performing Laboratory Location**

Mayo Clinic Laboratories - Rochester Main Campus

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.

# Test Definition: DMAID

Mycobacterium abscessus Complex Subspecies  
Identification, Next-Generation Sequencing

- 
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
  - Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

## Test Classification

This test was developed and its performance characteristics determined by Mayo Clinic in a manner consistent with CLIA requirements. It has not been cleared or approved by the US Food and Drug Administration.

## CPT Code Information

87153

## LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
DMAID	M abscessus subspecies ID, NGS	112278-7

Result ID	Test Result Name	Result LOINC® Value
DMAID	M abscessus subspecies ID, NGS	112278-7