

Overview

Method Name

High Performance Liquid Chromatography with Ultraviolet Detection; (HPLC-UV)

NY State Available

Yes

Specimen

Specimen Type

Varies

Specimen Required

Submit only 1 of the following specimens:

Plasma

Draw blood in a green-top sodium heparin tube(s), plasma gel tube is not acceptable. Spin down and send 3mL of plasma refrigerated in a plastic vial.

Serum

Draw blood in a plain, red-top tube(s), serum gel tube is not acceptable. Spin down and send 3 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.6 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	7 days	
	Ambient	72 hours	
	Frozen	180 days	

Clinical & Interpretive

Reference Values

Units: ng/mL

Expected fluvoxamine concentrations on recommended daily dosage regimens:  
50-900 ng/mL

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

3 to 7 days

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

80332

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FFVOX	Fluvoxamine (Luvox)	10988-4

Result ID	Test Result Name	Result LOINC® Value
FFVOX	Fluvoxamine (Luvox)	10988-4