

Overview

Method Name

Enzymatic method

NY State Available

No

Specimen

Specimen Type

Body Fluid

Specimen Required

Specimen Type: Tracheal or bronchial fluid

Sources: Tracheal or bronchial fluid

Container/Tube: Standard transport tube

Specimen Volume: 1 mL

Collection Instructions:

- 1. Collect 1 mL tracheal or bronchial fluid.
- 2. Transfer to a tightly capped plastic tube and freeze immediately.
- 3. Ship frozen.

Specimen Minimum Volume

0.5 mL

Reject Due To

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Body Fluid	Frozen	30 days	

Clinical & Interpretive

Reference Values

Protein Unit: mg/mL

Pepsin A Unit: ng/mL

Pepsin A Reference Range (ng/mL):

<12.5 negative

12.5-100 weak to moderate positive

>100 strong positive

Performance

PDF Report

No

Day(s) Performed

Monday through Thursday

Report Available

4 to 9 days

Performing Laboratory Location

Nemours Children's Health-Gastroenterology

Fees & Codes

- Fees
- Authorized users can sign in to [Test Prices](#) for detailed fee information.
 - Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
 - Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

83986

84157

82657

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FPEPA	Pepsin A Assay	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z3670	Specimen Type	31208-2
Z3671	pH	Not Provided
Z3672	Protein	Not Provided
Z3673	Pepsin A	Not Provided
Z3674	Interpretation	Not Provided