

**Overview****Method Name**

ImmunoCAP FEIA

**NY State Available**

Yes

**Specimen****Specimen Type**

Serum

**Specimen Required****Collection Container/Tube:** Red top**Acceptable:** Serum gel**Submission Container/Tube:** Plastic vial**Specimen Volume:** 0.5 mL of serum**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

|               |        |
|---------------|--------|
| Gross lipemia | Reject |
|---------------|--------|

**Specimen Stability Information**

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Serum         | Refrigerated (preferred) | 28 days  |                   |
|               | Frozen                   | 365 days |                   |
|               | Ambient                  | 28 days  |                   |

**Clinical & Interpretive****Reference Values**

&lt;0.35 kU/L

**Interpretation**

| <u>Class</u> | <u>IgE (kU/L)</u> | <u>Comment</u>       |
|--------------|-------------------|----------------------|
| 0            | <0.10             | Negative             |
| 0/1          | 0.10-0.34         | Equivocal/Borderline |
| 1            | 0.35-0.69         | Low Positive         |
| 2            | 0.70-3.49         | Moderate Positive    |
| 3            | 3.50-17.49        | High Positive        |
| 4            | 17.50-49.99       | Very High Positive   |
| 5            | 50.00-99.99       | Very High Positive   |
| 6            | >99.99            | Very High Positive   |

**Performance****PDF Report**

No

**Day(s) Performed**

Monday through Friday

**Report Available**

3 to 6 days

**Performing Laboratory Location**

Eurofins Viracor

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

[This test was developed, and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.](#)

**CPT Code Information**

86003

**LOINC® Information**

| Test ID | Test Order Name                 | Order LOINC® Value |
|---------|---------------------------------|--------------------|
| FSTAB   | Stachybotrys chartarum/atra IgE | 29616-0            |

| Result ID | Test Result Name                | Result LOINC® Value |
|-----------|---------------------------------|---------------------|
| Z4145     | Stachybotrys chartarum/atra IgE | 29616-0             |
| Z5805     | CLASS                           | Not Provided        |