

**Overview****Method Name**

Enzyme Immunoassay (EIA)

**NY State Available**

Yes

**Specimen****Specimen Type**

Serum

**Specimen Required****Specimen Type:** Serum**Collection Container/Tube:****Preferred:** Red**Acceptable:** Serum gel**Submission Container/Tube:** Plastic vial**Specimen Volume:** 0.5 mL Serum**Collection Instructions:**

1. Centrifuge and aliquot 0.5 mL of serum into a plastic vial.
2. Send refrigerated.

**Specimen Minimum Volume**

See Specimen Required

**Reject Due To**

|         |        |
|---------|--------|
| Lipemia | Reject |
|---------|--------|

**Specimen Stability Information**

| Specimen Type | Temperature              | Time     | Special Container |
|---------------|--------------------------|----------|-------------------|
| Serum         | Refrigerated (preferred) | 28 days  |                   |
|               | Ambient                  | 28 days  |                   |
|               | Frozen                   | 365 days |                   |

**Clinical & Interpretive**

**Clinical Information**

Refer to [www.eurofins-viracor.com/test-menu/](http://www.eurofins-viracor.com/test-menu/)

**Reference Values**

<0.35 kU/L

**Interpretation**

| Class | IgE (kU/L)  | Comment              |
|-------|-------------|----------------------|
| 0     | <0.35       | Below Detection      |
| 1     | 0.35-0.69   | Low Positive         |
| 2     | 0.70-3.49   | Moderate Positive    |
| 3     | 3.50-17.49  | Positive             |
| 4     | 17.50-49.99 | Strong Positive      |
| 5     | >49.99      | Very Strong Positive |

**Performance****PDF Report**

No

**Day(s) Performed**

Monday through Friday

**Report Available**

4 to 7 days

**Performing Laboratory Location**

Eurofins Viracor

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

[This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.](#)

**CPT Code Information**

---

86003

**LOINC® Information**

| Test ID | Test Order Name | Order LOINC® Value |
|---------|-----------------|--------------------|
| FRAD    | Radish IgE      | 7643-0             |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------|---------------------|
| Z3984     | Radish IgE       | 7643-0              |
| Z3985     | CLASS            | 102555-0            |