

## Overview

**Method Name**

Enzyme Immunoassay (EIA)

**NY State Available**

Yes

## Specimen

**Specimen Type**

Serum

**Specimen Required****Collection Container/Tube:****Preferred:** Red top**Acceptable:** Serum gel**Submission Container/Tube:** Plastic vial**Specimen Volume:** 0.5 mL Serum**Collection Instructions:**

1. Centrifuge and aliquot 0.5 mL of serum into a plastic vial.
2. Send refrigerate

**Specimen Minimum Volume**

See Specimen Required

**Reject Due To**

Lipemia	Reject
---------	--------

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

## Clinical & Interpretive

**Clinical Information**

Refer to [www.eurofins-viracor.com/test-menu](http://www.eurofins-viracor.com/test-menu)

**Reference Values**

<0.35 kU/L

**Interpretation**

Class	IgE (kU/L)	Comment
0	<0.35	Below Detection
1	0.35-0.69	Low Positive
2	0.70-3.49	Moderate Positive
3	3.50-17.49	Positive
4	17.50-49.99	Strong Positive
5	>49.99	Very Strong Positive

**Performance****PDF Report**

No

**Day(s) Performed**

Monday through Friday

**Report Available**

4 to 7 days

**Performing Laboratory Location**

Eurofins Viracor

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

**CPT Code Information**

86003

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
FASPE	Aspen (Populus tremuloides) IgE	7101-9

Result ID	Test Result Name	Result LOINC® Value
Z4012	Aspen (Populus tremuloides) IgE	7101-9
Z4013	CLASS	102152-6