

Test Definition: FRCE

Cedar Red (Juniperus virginiana) IgE

Overview

Method Name

ImmunoCAP FEIA

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Specimen Type: Serum
Container/Tube: Red or SST
Specimen Volume: 0.5 mL

Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is acceptable. Spin down and send 0.5 mL

of serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.5 mL

Reject Due To

Lipemia	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

Clinical & Interpretive

Reference Values

<0.35 kU/L



Test Definition: FRCE

Cedar Red (Juniperus virginiana) IgE

Interpretation

<u>Class</u>	IgE (kU/L)	<u>Comment</u>
0	<0.10	Negative
0/1	0.10-0.34	Equivocal/Borderline
1	0.35-0.69	Low Positive
2	0.70-3.49	Moderate Positive
3	3.50-17.49	High Positive
4	17.50-49.99	Very High Positive
5	50.0-99.99	Very High Positive
6	>99.99	Very High Positive

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

3 to 6 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

86003

LOINC® Information



Test Definition: FRCE

Cedar Red (Juniperus virginiana) IgE

Test ID	Test Order Name	Order LOINC® Value
FRCE	Cedar Red IgE	7195-1

Result ID	Test Result Name	Result LOINC® Value
Z4022	Cedar Red IgE	7195-1
Z4023	CLASS	102216-9