

Test Definition: FWWE

Hemp Western Water (Acnida tamariscina) IgE

Overview

Method Name

Enzyme Immunoassay (EIA)

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Collection Container/Tube:

Preferred: Red top **Acceptable**: Serum gel

Submission Container/tube: Plastic vial **Specimen Volume**: 0.5 mL serum

Collection Instructions:

1. Centrifuge and aliquot 0.5 mL of serum into a plastic vial.

2. Send refrigerate.

Specimen Minimum Volume

See Specimen Required

Reject Due To

Lipemia	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	28 days	
	Frozen	365 days	

Clinical & Interpretive

Clinical Information



Test Definition: FWWE

Hemp Western Water (Acnida tamariscina) IgE

https://www.eurofins-viracor.com

Reference Values

<0.35 kU/L

Interpretation

Class	IgE (kU/L)	Comment
0	<0.35	Below Detection
1	0.35-0.69	Low Positive
2	0.70-3.49	Moderate Positive
3	3.50-17.49	Positive
4	17.50-49.99	Strong Positive
5	>49.99	Very Strong Positive

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

4 to 7 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

86003



Test Definition: FWWE

Hemp Western Water (Acnida tamariscina) IgE

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FWWE	Hemp Western Water IgE	7400-5

Result ID Test Result Name		Result LOINC® Value	
Z4030	Hemp Western Water IgE	7400-5	
Z4031	CLASS	102366-2	