

## Overview

### Method Name

Only orderable as part of a profile. For more information see:

- MSP3 / Multiple Sclerosis (MS) Profile, Serum and Spinal Fluid
- OLIG / Oligoclonal Banding, Serum and Spinal Fluid

Collection Only

### NY State Available

No

## Specimen

### Specimen Type

Serum

### Specimen Required

Only orderable as part of a profile. For more information see:

- MSP3 / Multiple Sclerosis (MS) Profile, Serum and Spinal Fluid
- OLIG / Oligoclonal Banding, Serum and Spinal Fluid

**Specimen Type:** Serum

**Collection Container/Tube:**

**Preferred:** Serum gel

**Acceptable:** Red top

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

**Collection Instructions:**

1. Centrifuge and aliquot serum in plastic vial within 2 hours of collection.
2. Label specimen as serum.

### Specimen Minimum Volume

0.5 mL

### Reject Due To

Gross hemolysis	OK
Gross lipemia	OK

Gross icterus	OK
---------------	----

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	14 days	
	Ambient	14 days	
	Refrigerated	14 days	

**Clinical & Interpretive****Performance****PDF Report**

No

**Specimen Retention Time**

14 days

**Performing Laboratory Location**

Rochester

**Fees & Codes****Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

**Test Classification**

Not Applicable

**LOINC® Information**

Test ID	Test Order Name	Order LOINC® Value
XSRM	Additional sample for reflex OLIGS	No LOINC Needed

Result ID	Test Result Name	Result LOINC® Value
XSRM	Additional sample for Reflex OLIGS	No LOINC Needed