

**Overview****Method Name**

Only orderable as part of a profile. For more information see MSP3 / Multiple Sclerosis (MS) Cascade, Serum and Spinal Fluid.

Collection Only

**NY State Available**

No

**Specimen****Specimen Type**

Serum

**Specimen Required**

Only orderable as part of a profile. For more information see MSP3 / Multiple Sclerosis (MS) Cascade, Serum and Spinal Fluid.

**Collection Container/Tube:**

**Preferred:** Serum gel

**Acceptable:** Red top

**Submission Container/Tube:** Plastic vial

**Specimen Volume:** 1 mL

**Collection Instructions:**

1. Within 2 hours of collection, centrifuge and aliquot serum in plastic vial.
2. Label specimen as serum.

**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

Gross hemolysis	OK
Gross lipemia	OK
Gross icterus	OK

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	14 days	
	Ambient	14 days	
	Refrigerated	14 days	

## Clinical & Interpretive

## Performance

### PDF Report

No

### Specimen Retention Time

14 days

### Performing Laboratory Location

Mayo Clinic Laboratories - Rochester Superior Drive

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

Not Applicable

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
XSRM	Additional sample for reflex OLIGS	No LOINC Needed

Result ID	Test Result Name	Result LOINC® Value
XSRM	Additional sample for Reflex OLIGS	No LOINC Needed