

Overview

Testing Algorithm

This test helps to generate additional volume of whole blood to help complete testing in the scenario additional testing is needed.

Method Name

Collection Only

NY State Available

Yes

Specimen

Specimen Type

Whole blood

Specimen Required

Specimen Type: Whole blood

Container/Tube:

Preferred: Lavender top (EDTA)

Acceptable: Yellow top (ACD) or light-blue top (3.2% sodium citrate)

Specimen Volume: 3 mL

Collection Instructions:

1. Invert several times to mix blood.
2. Send specimen in original tube. **Do not** aliquot.

Specimen Stability Information

| Specimen Type | Temperature | Time | Special Container |
|---------------|---------------------|---------|-------------------|
| Whole blood | Ambient (preferred) | 14 days | |
| | Refrigerated | 14 days | |
| | Frozen | 14 days | |

Clinical & Interpretive

Performance

PDF Report

No

Performing Laboratory Location

Mayo Clinic Laboratories - Rochester Main Campus

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

Not Applicable

LOINC® Information

| Test ID | Test Order Name | Order LOINC® Value |
|---------|-----------------|--------------------|
| SC018 | Whole Blood | No LOINC Needed |

| Result ID | Test Result Name | Result LOINC® Value |
|-----------|------------------|---------------------|
| SC018 | Whole Blood | No LOINC Needed |