

SWI/SNF Chromatin Remodeling Complex Subunit ATPase 2 (SMARCA2 or BRM) Immunostain, Technical Component Only

# Overview

#### **Useful For**

Identifying and further subtyping SWI/SNF complex subunit loss associated with malignant rhabdoid tumors, epithelioid sarcomas, BRG1 (SMARCA4) deficient tumors, and other poorly differentiated malignancies

#### **Reflex Tests**

Test Id	Reporting Name	Available Separately	Always Performed
IHTOI	IHC Initial, Tech Only	No	No
IHTOA	IHC Additional, Tech Only	No	No

# **Testing Algorithm**

For the initial technical component only immunohistochemical (IHC) stain performed, the appropriate bill-only test ID will be reflexed and charged (IHTOI). For each additional technical component only IHC stain performed, an additional bill-only test ID will be reflexed and charged (IHTOA).

# **Method Name**

Immunohistochemistry (IHC)

#### **NY State Available**

Yes

# **Specimen**

#### **Specimen Type**

**TECHONLY** 

# **Ordering Guidance**

This test includes only technical performance of the stain (no pathologist interpretation is performed). If diagnostic consultation by a pathologist is required order PATHC / Pathology Consultation.

# **Shipping Instructions**

Attach the green "Attention Pathology" address label (T498) and the pink Immunostain Technical Only label included in the kit to the outside of the transport container.

# Specimen Required

Specimen Type: Tissue

Supplies: Immunostain Technical Only Envelope (T693)



SWI/SNF Chromatin Remodeling Complex Subunit ATPase 2 (SMARCA2 or BRM) Immunostain, Technical Component Only

Container/Tube: Immunostain Technical Only Envelope

**Submit:** 

Formalin-fixed, paraffin-embedded tissue block

OR

2 Unstained, positively charged glass slides (25- x 75- x 1-mm) per test ordered; sections 4-microns thick

# **Digital Image Access**

- 1. Information on accessing digital images of immunohistochemical (IHC) stains and the manual requisition form can be accessed through this website: <a href="https://news.mayocliniclabs.com/pathology/digital-imaging/">https://news.mayocliniclabs.com/pathology/digital-imaging/</a>
- 2. Clients ordering stains using a manual requisition form will not have access to digital images.
- 3. Clients wishing to access digital images must place the order for IHC stains electronically. Information regarding digital imaging can be accessed through this website: <a href="https://news.mayocliniclabs.com/pathology/digital-imaging/#section3">https://news.mayocliniclabs.com/pathology/digital-imaging/#section3</a>

#### **Forms**

If not ordering electronically, complete, print, and send a <u>Immunohistochemical (IHC)/In Situ Hybridization (ISH) Stains</u>
Request (T763) with the specimen.

# Reject Due To

Wet/frozen	Reject
tissue	
Cytology	Reject
smears	
Nonformalin	Reject
fixed tissue	
Nonparaffin	Reject
embedded	
tissue	
Noncharged	Reject
slides	
ProbeOn slides	Reject
Snowcoat	Reject
slides	

# **Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
TECHONLY	Ambient (preferred)		
	Refrigerated		

# **Clinical & Interpretive**



SWI/SNF Chromatin Remodeling Complex Subunit ATPase 2 (SMARCA2 or BRM) Immunostain, Technical Component Only

#### **Clinical Information**

BRM is one of two mutually exclusive catalytic ATPase subunits of the SWI/SNF complex. The loss of at least one SWI/SNF subunit has been associated with undifferentiated carcinomas and/or advanced tumorigenesis. The homogenous and heterogenous loss of BRM can occur independently or concordantly with other SWI/SNF subunits.

# Interpretation

This test does not include pathologist interpretation, only technical performance of the stain. If interpretation is required, order PATHC / Pathology Consultation for a full diagnostic evaluation or second opinion of the case.

The positive and negative controls are verified as showing appropriate immunoreactivity.

Interpretation of this test should be performed in the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

#### **Cautions**

Age of a cut paraffin section can affect immunoreactivity. Stability thresholds vary widely among published literature and are antigen dependent. Best practice is for paraffin sections to be cut within 6 weeks.

The charge of glass slides can be affected by environmental factors and subsequently may alter slide staining. Sending unsuitable glass slides can result in inconsistent staining due to poor slide surface chemistry.

Best practices for storage of positively charged slides:

- -Minimize time slides are stored after being unpackaged
- -Limit exposure to high humidity and heat
- -Minimize exposure to plastics

#### **Clinical Reference**

- 1. Marquez-Vilendrer SB, Rai SK, Gramling SJ, Lu L, Reisman DN. Loss of the SWI/SNF ATPase subunits BRM and BRG1 drives lung cancer development. Oncoscience. 2016;3(11-12):322-336
- 2. Jancewicz I, Siedlecki1 JA, Sarnowski TJ, Sarnowska E. BRM: the core ATPase subunit of SWI/SNF chromatin-remodelling complex-a tumour suppressor or tumour-promoting factor?. Epigenetics & Chromatin. 2019;12(1):68
- 3. Chang B, Sheng W, Wang L, et al. SWI/SNF Complex-deficient Undifferentiated Carcinoma of the Gastrointestinal Tract: Clinicopathologic Study of 30 Cases With an Emphasis on Variable Morphology, Immune Features, and the Prognostic Significance of Different SMARCA4 and SMARCA2 Subunit Deficiencies. Am J Surg Pathol. 2022;46(7):889-906 4. Li L, Fan Xs, Xia Qy, et al. Concurrent loss of INI1, PBRM1, and BRM expression in epithelioid sarcoma: implications for the cocontributions of multiple SWI/SNF complex members to pathogenesis. Hum Pathol. 2014;45(11):2247-2254

#### **Performance**

# **Method Description**



SWI/SNF Chromatin Remodeling Complex Subunit ATPase 2 (SMARCA2 or BRM) Immunostain, Technical Component Only

Immunohistochemistry on sections of paraffin-embedded tissue.(Unpublished Mayo method)

# PDF Report

No

# Day(s) Performed

Monday through Friday

## **Report Available**

1 to 3 days

# **Specimen Retention Time**

Until staining is completed

#### **Performing Laboratory Location**

Mayo Clinic Laboratories - Rochester Main Campus

# **Fees & Codes**

### Fees

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact Customer Service.

#### **Test Classification**

This test has been cleared, approved, or is exempt by the US Food and Drug Administration and is used per manufacturer's instructions. Performance characteristics were verified by Mayo Clinic in a manner consistent with CLIA requirements.

# **CPT Code Information**

88342-TC, Primary 88341-TC, If additional IHC

#### **LOINC®** Information

Test ID	Test Order Name	Order LOINC® Value
HBRM	BRM (SMARCA2) IHC, Tech Only	No LOINC Needed

Result ID	Test Result Name	Result LOINC® Value
622928	BRM (SMARCA2) IHC, Tech Only	No LOINC Needed