

# **Test Definition: FMARC**

Bupivacaine (Marcaine)

# **Overview**

### **Method Name**

Liquid Chromatography/ Tandem Mass Spectrometry (LC/MS/MS)

#### **NY State Available**

Yes

# **Specimen**

# **Specimen Type**

Serum Red

## **Specimen Required**

Draw blood in a plain red-top tube(s), **serum gel tube is not acceptable**. Spin down and send 1 mL of serum frozen in a plastic vial.

## **Specimen Minimum Volume**

0.5 mL

### **Reject Due To**

Hemolysis	NA NA
Lipemia	NA
Icterus	NA
Other	NA

## **Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum Red	Frozen (preferred)	180 days	
	Ambient	72 hours	
	Refrigerated	72 hours	

# **Clinical & Interpretive**

# **Reference Values**

Reference Range: 0.1 - 4.0 ug/mL



# **Test Definition: FMARC**

Bupivacaine (Marcaine)

## **Performance**

### **PDF Report**

No

## Day(s) Performed

Monday through Sunday

### **Report Available**

5 to 11 days

## **Performing Laboratory Location**

Medtox Laboratories, Inc.

### **Fees & Codes**

### **Fees**

- Authorized users can sign in to <u>Test Prices</u> for detailed fee information.
- Clients without access to Test Prices can contact <u>Customer Service</u> 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact <u>Customer Service</u>.

# **Test Classification**

This test was developed and its performance characteristics determined by LabCorp. It has not been cleared or approved by the Food and Drug Administration.

## **CPT Code Information**

80299

# **LOINC®** Information

Test ID	Test Order Name	Order LOINC® Value
FMARC	Bupivacaine (Marcaine)	3412-4

Result ID	Test Result Name	Result LOINC® Value
Z5052	Bupivacaine	3412-4