

**Overview****Method Name**

Enzyme Immunoassay (EIA)

**NY State Available**

Yes

**Specimen****Specimen Type**

Serum

**Specimen Required****Specimen Type:** Serum**Container/Tube:** Red or SST**Specimen Volume:** 0.5 mL**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 0.5 mL of serum refrigerated in a plastic vial.**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

Lipemia	Reject
---------	--------

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Ambient	28 days	
	Refrigerated (preferred)	28 days	
	Frozen	365 days	

**Clinical & Interpretive****Reference Values**

&lt;0.35 kU/L

**Interpretation**

---

<u>Class</u>	<u>IgE (kU/L)</u>	<u>Comment</u>
0	<0.35	Below Detection
1	0.35-0.69	Low Positive
2	0.70-3.49	Moderate Positive
3	3.50-17.49	Positive
4	17.50-49.99	Strong Positive
5	>49.99	Very Strong Positive

## Performance

### Method Description

#### PDF Report

No

#### Day(s) Performed

Monday through Friday

#### Report Available

4 to 7 days

#### Performing Laboratory Location

Eurofins Viracor

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

### CPT Code Information

86003

### LOINC® Information

---

Test ID	Test Order Name	Order LOINC® Value
FPINP	Pine Ponderosa IgE	7610-9

Result ID	Test Result Name	Result LOINC® Value
Z5335	Pine Ponderosa IgE	7610-9
Z5336	CLASS	Not Provided