

Overview

Method Name
Reversed-phase UPLC-MS/MS

NY State Available
Yes

Specimen

Specimen Type
Serum

Specimen Required
Specimen Type: Serum
Container/Tube: Red
Specimen Volume: 1 mL
Collection Instructions: Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Allow to clot for 30 to 60 minutes, spin down and send 1 mL of serum Frozen in a sterile, screw top tube.

Specimen Minimum Volume
0.5 mL

Reject Due To

Gross hemolysis:	Reject
Thawing:	Warm Reject; Cold OK
Gross lipemia:	Reject
Other:	SST

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	30 days	
	Ambient	14 days	
	Refrigerated	14 days	

Clinical & Interpretive

Reference Values

The analytical measuring range is 0.1-10 mcg/mL. The therapeutic range has not yet been established.

Performance

PDF Report

No

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

80299

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FIVCZ	Isavuconazole	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FIVCZ	Isavuconazole	Not Provided