

## Overview

**Method Name**

Reversed-phase UPLC-MS/MS

**NY State Available**

Yes

## Specimen

**Specimen Type**

Serum

**Specimen Required****Specimen Type:** Serum**Container/Tube:** Red**Specimen Volume:** 1 mL**Collection Instructions:** Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Allow to clot for 30 to 60 minutes, spin down and send 1 mL of serum Frozen in a sterile, screw top tube.**Specimen Minimum Volume**

0.5 mL

**Reject Due To**

Gross hemolysis:	Reject
Thawing:	Warm Reject; Cold OK
Gross lipemia:	Reject
Other:	SST

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Frozen (preferred)	30 days	
	Ambient	14 days	
	Refrigerated	14 days	

## Clinical & Interpretive

## Reference Values

The analytical measuring range is 0.1-10 mcg/mL. The therapeutic range has not yet been established.

## Performance

### PDF Report

No

### Performing Laboratory Location

Eurofins Viracor

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### Test Classification

This test was developed and its performance characteristics determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

### CPT Code Information

80299

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FIVCZ	Isavuconazole	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FIVCZ	Isavuconazole	Not Provided