

Overview

Method Name

Quantitative Spectrophotometry

NY State Available

Yes

Specimen

Specimen Type

CSF

Specimen Required

Collection Container/Tube: Leak-proof container

Submission Container/Tube: Sterile vial

Specimen Volume: 0.5 mL

Collection Instructions:

1. Centrifuge cerebrospinal fluid specimen at ambient temperature.
2. Transfer 0.5 mL fluid (supernatant) into a sterile vial and freeze.
3. Send frozen. **Specimen must remain frozen** until received at performing lab.
4. Indicate source on requisition.

Specimen Minimum Volume

0.2 mL

Reject Due To

Bronchoalveolar lavage (BAL)	Reject
Glass tubes	Reject
Turbid specimens	Reject
Whole blood	Reject

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
CSF	Frozen (preferred)	30 days	
	Ambient	24 hours	
	Refrigerated	7 days	

Clinical & Interpretive**Clinical Information**

Refer to <https://ltd.aruplab.com/>

Reference Values

0-9 U/L

Performance**PDF Report**

No

Day(s) Performed

Tuesday, Thursday, Sunday

Report Available

1 to 8 days

Performing Laboratory Location

ARUP Laboratories

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by ARUP Laboratories. It has not been cleared or approved by the US Food and Drug Administration. This test was performed in a CLIA certified laboratory and is intended for clinical purposes.

CPT Code Information

84311

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
---------	-----------------	--------------------

FADCF	Adenosine Deaminase in CSF	35703-8
Result ID	Test Result Name	Result LOINC® Value
Z5892	Adenosine Deaminase, CSF	35703-8