

Overview

Special Instructions

- [Versiti Platelet Neutrophil Immunology Lab Form](#)

Highlights

This assay includes: Testing for antibodies against HPA-1a/b, HPA-2a/b, HPA-3a/b, HPA-4a, HPA-5a/b, GPIIb/IIIa, GPIa/IIa, GPIb/IX. GPIV and Class I HLA.

Method Name

Platelet Antibody Bead Array (PABA) and Flow Cytometry

NY State Available

Yes

Specimen

Specimen Type

Serum

Specimen Required

Specimen Type: Serum

Collection Container/Tube: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 5 mL

Collection Instruction:

- 1.Centrifuge and aliquot 5 mL serum into a plastic vial.
- 2.Send refrigerated.

Forms

Complete Platelet Neutrophil Immunology Laboratory form

Specimen Minimum Volume

3 mL

Reject Due To

Gross hemolysis	Reject
-----------------	--------

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	7 days	
	Frozen	365 days	

Clinical & Interpretive

Reference Values

Negative

Performance

PDF Report

Referral

Day(s) Performed

Varies

Report Available

9 to 11 days

Performing Laboratory Location

Versiti Wisconsin, Inc.

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Versiti Wisconsin, Inc. It has not been cleared or approved by the US Food and Drug Administration. This test is used for clinical purposes. It should not be regarded as investigational or for research. This laboratory is certified under the Clinical Laboratory Improvement Amendments (CLIA) as qualified to perform high complexity clinical laboratory testing.

CPT Code Information

86022

Test Definition: FNNAT

Neonatal Alloimmune Thrombocytopenia
(NAIT)-Serial Monitoring Maternal only

86023 x 2

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FNNAT	Neonatal Alloimmune Thromb (NAIT)	Not Provided

Result ID	Test Result Name	Result LOINC® Value
FNNAT	Neonatal Alloimmune Thromb (NAIT)	Not Provided