

## Overview

**Method Name**

Fluoroimmunoassay Multiplex Bead Assay

**NY State Available**

Yes

## Specimen

**Specimen Type**

Serum

**Specimen Required****Collection Container/Tube:****Preferred:** Red top**Acceptable:** Serum gel**Submission Container/tube:** Plastic vial**Specimen Volume:** 1 mL Serum**Collection Instructions:**

1. Allow blood to clot for 30 to 60 minutes.
2. Centrifuge at ambient temperature, 1000 x g for 15 minutes.
3. Aliquot 1 mL of serum into a plastic vial.
4. Send refrigerated.

**Specimen Minimum Volume**

Serum: 0.1 mL

**Reject Due To**

All specimens will be evaluated by the processing and performing laboratories for test suitability

**Specimen Stability Information**

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	42 days	
	Ambient	42 days	
	Frozen		

---

**Clinical & Interpretive****Clinical Information**

Refer to [www.eurofins-viracor.com/test-menu/](http://www.eurofins-viracor.com/test-menu/)

**Reference Values**

Pneumococcal Antibody Panel - PCV15

Pneumo Ab Type 1: >1.3 ug/mL

Pneumo Ab Type 3: >1.3 ug/mL

Pneumo Ab Type 4: >1.3 ug/mL

Pneumo Ab Type 5: >1.3 ug/mL

Pneumo Ab Type 6A: >1.3 ug/mL

Pneumo Ab Type 6B: >1.3 ug/mL

Pneumo Ab Type 7F: >1.3 ug/mL

Pneumo Ab Type 9V: >1.3 ug/mL

Pneumo Ab Type 14: >1.3 ug/mL

Pneumo Ab Type 18C: >1.3 ug/mL

Pneumo Ab Type 19A: >1.3 ug/mL

Pneumo Ab Type 19F: >1.3 ug/mL

Pneumo Ab Type 22F: >1.3 ug/mL

Pneumo Ab Type 23F: >1.3 ug/mL

Pneumo Ab Type 33F: >1.3 ug/mL

**Performance****PDF Report**

No

**Day(s) Performed**

Monday through Friday

**Report Available**

9 to 11 days

**Performing Laboratory Location**

Eurofins Viracor

**Fees & Codes**

## Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

## Test Classification

This test was developed and its performance characteristics were determined by Eurofins Viracor. It has not been cleared or approved by the U.S. Food and Drug Administration.

## CPT Code Information

86581

## LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FPA15	Pneumococcal Antibody Panel - PCV15	Not Provided

Result ID	Test Result Name	Result LOINC® Value
Z6409	Pneumo Ab Type 1	85955-3
Z6410	Pneumo Ab Type 3	86081-7
Z6411	Pneumo Ab Type 4	86108-8
Z6412	Pneumo Ab Type 5	86129-4
Z6413	Pneumo Ab Type 6A	86138-5
Z6414	Pneumo Ab Type 6B	40905-2
Z6415	Pneumo Ab Type 7F	40911-0
Z6416	Pneumo Ab Type 9V	40926-8
Z6417	Pneumo Ab Type 14	85992-6
Z6418	Pneumo Ab Type 18C	40913-6
Z6419	Pneumo Ab Type 19A	40915-1
Z6420	Pneumo Ab Type 19F	86021-3
Z6421	Pneumo Ab Type 22F	86053-6
Z6422	Pneumo Ab Type 23F	86061-9
Z6423	Pneumo Ab Type 33F	47349-6