

Overview

Method Name

Liquid Chromatography/Tandem Mass Spectrometry (LC/MS/MS)

NY State Available

Yes

Specimen

Specimen Type

Serum Red

Specimen Required

Collection Container/Tube: Red top

Submission Container/Tube: Plastic vial

Specimen Volume: 3 mL

Collection Instructions: Draw blood in a plain, red-top tube(s), **serum-gel tube is not acceptable**. Centrifuge and send 3 mL of serum refrigerate in a plastic vial.

Specimen Minimum Volume

0.5 mL

Reject Due To

Hemolysis	Mild reject; Gross reject
Lipemia	Mild reject; Gross reject
Icterus	Mild reject; Gross reject
Other	Specimens other than Serum collected in red top. Plasma or whole blood collected in Sodium Heparin or EDTA. Serum gel tubes or Plasma gel are NOT acceptable.

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum Red	Ambient	72 hours	
	Refrigerated (preferred)	14 days	
	Frozen	180 days	

Clinical & Interpretive

Reference Values

Reference Range: 1.0 - 10 ng/mL

Performance**PDF Report**

No

Day(s) Performed

Monday through Sunday

Report Available

5 to 9 days

Specimen Retention Time

2 weeks

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees & Codes**Fees**

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

80342

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
PROLX	Fluphenazine (Prolixin)	3650-9

Result ID	Test Result Name	Result LOINC® Value
Z1122	Fluphenazine	3650-9