

Overview

Method Name
High Performance Liquid Chromatography with Fluorescence Detection (HPLC-FL)

NY State Available
Yes

Specimen

Specimen Type
Varies

Specimen Required
Submit only 1 of the following specimens

Plasma
Draw blood in a green-top (sodium heparin) tube(s), plasma gel tube is not acceptable. Spin down and send 1 mL sodium heparin plasma refrigerated in a plastic vial.

Serum
Draw blood in a plain red-top tube(s), serum gel tube is not acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume
0.25 mL

Reject Due To

Hemolysis	NA
Lipemia	NA
Icterus	NA
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Ambient	72 hours	
	Frozen	180 days	

Clinical & Interpretive

Reference Values

Units: ug/mL

Therapeutic and toxic ranges have not been established.

Expected serum vigabatrin concentrations in patients receiving recommended daily dosages: 20 – 160 ug/mL

Performance

PDF Report

No

Day(s) Performed

Monday through Sunday

Report Available

2 to 6 days

Performing Laboratory Location

Medtox Laboratories, Inc.

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

CPT Code Information

80339

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
VIGA	Vigabatrin (Sabril)	30042-6

Result ID	Test Result Name	Result LOINC® Value
Z1176	Vigabatrin (Sabril)	30042-6