

## Overview

### Method Name

Gas Chromatography-Mass Spectrometry (GC/MS)

### NY State Available

Yes

## Specimen

### Specimen Type

Varies

### Specimen Required

Submit only 1 of the following specimens:

#### Plasma

Draw blood in a green-top (sodium heparin) tube(s), **plasma gel tube is not acceptable**. Spin down and send 5 mL sodium heparin plasma refrigerated in a plastic vial.

#### Serum

Draw blood in a plain red-top tube(s), **serum gel tube is not acceptable**. Spin down and send 5 mL of serum refrigerated in a plastic vial.

### Specimen Minimum Volume

2 mL

### Reject Due To

Hemolysis:	NA
Thawing:	Warm OK: Cold OK
Lipemia:	NA
Icterus:	NA
Other	NA

### Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Varies	Refrigerated (preferred)	14 days	
	Ambient	48 hours	
	Frozen	365 days	

## Clinical & Interpretive

### Reference Values

Units: Flunitrazepam      ng/mL  
7-Aminoflunitrazepam ng/mL

Peak plasma Flunitrazepam concentrations in patients receiving chronic, recommended dosages: 10-20 ng/mL.

Note: Flunitrazepam is not legally marketed in the United States.

## Performance

### PDF Report

No

### Day(s) Performed

Monday through Sunday

### Report Available

5 to 9 days

### Performing Laboratory Location

Medtox Laboratories, Inc.

## Fees & Codes

### Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

### CPT Code Information

80346

### LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FFLRO	Flunitrazepam (Rohypnol)	Not Provided

Result ID	Test Result Name	Result LOINC® Value
-----------	------------------	---------------------

---

Z1075	Flunitrazepam	73828-6
Z1074	7-Aminoflunitrazepam	73829-4