

Overview

Method Name

ImmunoCAP

NY State Available

No

Specimen

Specimen Type

Serum

Specimen Required

Draw blood in a plain red-top tube(s), serum gel tube(s) is acceptable. Spin down and send 1 mL of serum refrigerated in a plastic vial.

Specimen Minimum Volume

0.5 mL

Reject Due To

Hemolysis	Mild OK; Gross Reject
Lipemia	Mild OK; Gross Reject
Icterus	Mild OK; Gross Reject
Other	NA

Specimen Stability Information

Specimen Type	Temperature	Time	Special Container
Serum	Refrigerated (preferred)	28 days	
	Ambient	7 days	
	Frozen	365 days	

Clinical & Interpretive

Reference Values

Reference ranges have not been established for food-specific IgG4 tests. The clinical utility of food-specific IgG4 tests has not been clearly established. These tests can be used in special clinical situations to select foods for evaluation by diet.

elimination and challenge in patients who have food-related complaints, and to evaluate food allergic patients prior to food challenges. The presence of food-specific IgG4 has been studied in response to various oral food immunotherapy treatments but cutoffs have not been established.

Performance

PDF Report

No

Day(s) Performed

Monday through Friday

Report Available

3 to 7 days

Performing Laboratory Location

Eurofins Viracor

Fees & Codes

Fees

- Authorized users can sign in to [Test Prices](#) for detailed fee information.
- Clients without access to Test Prices can contact [Customer Service](#) 24 hours a day, seven days a week.
- Prospective clients should contact their account representative. For assistance, contact [Customer Service](#).

Test Classification

This test was developed and its performance characteristics determined by Viracor Eurofins. It has not been cleared or approved by the U.S. Food and Drug Administration.

CPT Code Information

86001 x 8

LOINC® Information

Test ID	Test Order Name	Order LOINC® Value
FFPG4	Food Panel IgG4 (532)	Not Provided
Result ID	Test Result Name	Result LOINC® Value
Z4353	Corn IgG4	56231-4
Z4354	Egg White IgG4	56253-8
Z4355	Egg Yolk IgG4	56254-6

Z4356	Milk Cow IgG4	56234-8
Z4357	Peanut IgG4	56351-0
Z4358	Soybean IgG4	56390-8
Z4359	Wheat IgG4	56417-9
Z4360	Yeast(<i>Saccharomyces cerevisiae</i>)IgG4	56449-2